

The background of the slide features a blurred image of a laboratory microplate. A pipette is shown in the upper right corner, dispensing a drop of bright pink liquid into one of the wells. The overall scene is brightly lit, with a soft, out-of-focus background.

PacBio

HiFi sequencing and software v11.0 release: Technical overview for Sequel II and Sequel IIe system users

Sequel II and IIe Systems ICS v11.0 / SMRT Link v11.0

PN 102-399-900 Version 01 (April 2022)

HiFi sequencing and software v11.0 release

Technical overview

A. New consumables for SMRTbell library preparation & sequencing preparation (ABC)

- SMRTbell prep kit 3.0 (SPK 3.0)
- AMPure PB bead size selection kit
- Binding kit 3.1 & cleanup beads
- Binding kit 3.2 & cleanup beads

B. SPK 3.0 WGS sample preparation updates

- SPK 3.0 WGS workflow overview
- SPK 3.0 WGS example performance data

C. SMRT Link Sample Setup updates

- New High-throughput (HT) Sample Setup mode overview
- Sample setup HT mode annealing / binding / cleanup (ABC) and DNA internal control dilution procedure

D. SMRT link run design updates

- GUI updates and new default CCS analysis output changes
- New 5mC detection option
- New heteroduplex detection option
- New AAV application type option
- New on-instrument demultiplexing option

E. SMRT Link Run QC updates

- New instrument status view
- New 5mC-specific report plots

F. SMRT Link SMRT Analysis updates

- Updated GUI nomenclature
- Updated microbial genome analysis application

G. APPENDIX: Technical documentation & applications support resources

- Sequel IIe system documentation
- SMRT Link & other data analysis documentation
- Application technical overviews
- DNA extraction literature
- Sample preparation literature
- Example PacBio data sets



Summary overview of key features, improvements & customer benefits

System v11.0 key feature updates



ICS v11.0

Sequel II and IIe systems only



SMRT Link v11.0

Sequel II and IIe systems only

Workflow step	What's New	Notes
DNA sample prep	 <ul style="list-style-type: none"> Human whole blood DNA extraction tech note 	<ul style="list-style-type: none"> Additional DNA extraction guidance coming
SMRTbell library prep	  <ul style="list-style-type: none"> SMRTbell prep kit 3.0 / AMPure PB bead size selection kit Updated protocols (WGS / amplicon / HiFiViral / Iso-Seq method) and new AAV protocol 	<ul style="list-style-type: none"> New reagents to accelerate, unify & streamline library prep TPK 1.0 & TPK 2.0 still available
Sample setup (ABC)	  <ul style="list-style-type: none"> Sequel II binding kit 3.1 & 3.2 Sample Setup HT mode 	<ul style="list-style-type: none"> Uses Polymerases 2.1/2.2 in a new kit configuration to unify/streamline ABC Binding kits 2.1/2.2 still available
Run design	  <ul style="list-style-type: none"> On-instrument 5mC / heteroduplex detection (HD) / demultiplexing AAV sequencing mode 	<ul style="list-style-type: none"> 5mC, HD and demultiplexing capabilities are also available in SMRT Link for Sequel II system
Sequencing	  <ul style="list-style-type: none"> SMRT Cell 8M single-use tray Sequel II sequencing kit 2.0 (1 rxn) 	<ul style="list-style-type: none"> Functionally the same basic parts as before with updated kit configurations to enable increased run flexibility
Run QC	 <ul style="list-style-type: none"> Instrument status monitoring New 5mC report plots 	<ul style="list-style-type: none"> Provides real-time ZMW loading performance information during sequencing runs
SMRT Analysis	 <ul style="list-style-type: none"> Refreshed GUI emphasizes HiFi data Microbial genome analysis application 	<ul style="list-style-type: none"> Deprecates obsoleted subread (CLR) data-based workflows

System v11.0 release key benefits

New products enable simplification & acceleration of sample preparation workflows and consumables ordering

New consumables for SMRTbell library preparation & sequencing preparation


- **New streamlined reagent kit products** for SMRTbell library construction and annealing/binding/cleanup (ABC) of constructed libraries for HiFi sequencing
- Improved consolidation of required core template prep reagents into fewer product bundles reduces complexity and allows for **easier consumables ordering**

Updated SMRTbell library construction workflows to support SPK 3.0

- New whole genome sequencing (WGS) library prep protocol using SPK 3.0 enables **reduced workflow times** and **reduced DNA input** requirements; and **eliminates gel-based size selection requirements** for constructing large/small genome WGS and shotgun metagenomic libraries
- Amplicon, Iso-Seq method, and HiFiViral SARS-CoV-2 library prep protocols are updated to support SPK 3.0 and enable **simpler, more unified sample prep workflows** across different applications
- New adeno-associated virus (AAV) library prep protocol using SPK 3.0 support biopharma and **AAV gene editing vector R&D**

New single-reaction consumables products

- **New Sequel II sequencing kit 2.0 (1 rxn)** and **SMRT Cell 8M single use tray** products
- Enables **improved run schedule flexibility** for lower-throughput customers (no need to wait for 4 samples to accrue) or smaller project sizes



System v11.0 release key benefits (cont.)

New software enables better data, faster – and provides push-button access to the epigenome


Updated software for Sequel II and IIe systems

ICS v11.0



- ICS v11.0 is available for Sequel II and IIe systems only

SMRT Link v11.0



- SMRT Link v11.0 requires Sequel II and IIe systems to be upgraded to ICS v11.0

SMRT Link v11.0 Feature updates



Sample Setup

New High-Throughput mode feature provides a more streamlined workflow to efficiently process single samples or multiple samples in parallel using automation



Run Design

New on-instrument 5mC calling in CpG motifs.
New heteroduplex detection option detects and resolves heteroduplex reads
New on-instrument demultiplexing option
New adeno-associated virus (AAV) application type option



Run QC

New Instrument Status view in Run QC provides real-time ZMW loading performance information about PacBio instruments that are connected to SMRT Link




SMRT Analysis

New Microbial Genome Analysis application in SMRT Analysis for *de novo* assembly and base modification detection analysis of microbial genomes using HiFi reads

PacBio
SMRT® Link

General usability and user experience improvements to SMRT Link graphical user interface to emphasize HiFi sequencing, QC and analysis and deprecate obsolete subread (CLR) data-based workflows



New consumables for SMRTbell library prep and sequencing prep (ABC)

New consumables for SMRTbell library preparation


New products enable simplification & acceleration of sample preparation workflows and consumables ordering

New consumables for SMRTbell library preparation

New consumables for library construction			New consumables for size selection		
SMRTbell prep kit 3.0 Contains sufficient reagents to perform 24 SMRTbell library preparation reactions		102-182-700	AMPure PB beads size selection kit Contains sufficient reagents to perform 48 bead-based size selection reactions		102-182-500
	SMRTbell prep kit 3.0 <i>Includes all core reagents needed for SMRTbell library construction</i>			AMPure PB beads [5 mL] <i>Used for performing bead-based library size selection of WGS SMRTbell libraries</i>	
	SMRTbell cleanup beads [10 mL] <i>Used for routine DNA purification during SMRTbell library construction</i>			Elution buffer [50 mL] <i>Used for preparing diluted AMPure PB beads for size selection</i>	
	Low TE buffer [10 mL] <i>Used for shearing genomic DNA samples for WGS SMRTbell library construction</i>				

SMRTbell prep kit 3.0 configuration


- SPK 3.0 (102-182-700) contains all reagents needed to perform 24 SMRTbell library preparation reactions



Note: Primer for annealing is now included in Sequel II Binding Kit 3.1/3.2 to streamline ABC workflows.

AMPure PB beads size selection kit

- Bundle of **AMPure PB beads** and **Elution buffer**
- Depletes library fragments <5 kb
- Faster, more scalable, and lower capital equipment costs than automated gel-based size selection tools



Recommended








	Size-selection Beads	Gel Cassettes
Size Cutoff	Fixed (>5 kb)	✓ Adjustable
Hands-on Time	✓ Low	✓ Low
Run Time	✓ Minutes	Hours
Automation	✓ Standard liquid handler	Separate instrument
Cost	✓ Low	High
DNA Recovery	✓ High	Low
DNA Input to Library Prep per SMRT® Cell	✓ 1 µg	1.7 µg


Note: Gel size selection is still supported (adds 2 or more hours of run time and increases DNA input requirements). Refer to **Technical Note: Alternative Size Selection Methods for SMRTbell prep kit 3.0** ([102-326-503](#)), which provides detailed guidance for size selection of WGS libraries using automated gel-based DNA size selection tools.

New consumables for sequencing preparation (ABC)

New products enable simplification & acceleration of sequencing preparation workflows and consumables ordering

New consumables for sequencing preparation

New binding kit consumables		New SMRT Cell and sequencing kit consumables		
Binding kit 3.1 and cleanup beads Contains sufficient reagents to perform 24 binding reactions with Polymerase 3.1 for samples with inserts <3 kb		102-333-400	 SMRT Cell 8M single use tray Contains 1 SMRT Cell to be used with the Sequel II and IIe Systems.	102-182-500
	Sequel II binding kit 3.1 and DNA internal control 3.1 <i>Include Sequencing primer 3.1, Polymerase 2.1 and DNA internal control 3.1</i>		 Sequel II sequencing kit 2.0 (1 rxn) Contains sufficient reagents to support sequencing on one SMRT Cell 8M	102-194-400
	SMRTbell cleanup beads [10 mL] <i>Used for SMRTbell template cleanup after polymerase binding step</i>			
Binding kit 3.2 and cleanup beads Contains sufficient reagents to perform 24 binding reactions with Polymerase 3.2 for samples with inserts ≥3 kb		102-333-300		
	Sequel II binding Kit 3.2 and DNA internal control 3.2 <i>Include Sequencing primer 3.2, Polymerase 2.2 and DNA internal control 3.2</i>			
	SMRTbell cleanup beads [10 mL] <i>Used for SMRTbell template cleanup after polymerase binding step</i>			



Summary comparison of new System v11.0 release consumables versus previous consumables

New System v11.0 consumables provide an enhanced user experience with simplified, unified, scalable workflows

New v11.0 product	Number of reactions	Part number	Replaces*
SMRTbell prep kit 3.0	24	102-182-700	<ul style="list-style-type: none"> Express TPK 2.0 Enzyme cleanup kit 2.0 AMPure PB beads/ProNex
Binding kit 3.1 and cleanup beads	24 (up to 96 SMRT Cells 8M)	102-333-400	<ul style="list-style-type: none"> Binding kit 2.1 and internal control 1.0 AMPure PB beads/ProNex beads
Binding kit 3.2 and cleanup beads	24 (up to 96 SMRT Cells 8M)	102-333-300	<ul style="list-style-type: none"> Binding kit 2.2 and internal control 1.0 AMPure PB beads/ProNex beads Sequencing primer v5
AMPure PB beads size-selection kit	48	102-182-500	<ul style="list-style-type: none"> Gel-based size-selection
Sequel II sequencing kit 2.0 (1 rxn)	1	102-194-400	<ul style="list-style-type: none"> N/A
SMRT Cell 8M single-use tray	1	102-281-700	<ul style="list-style-type: none"> N/A

* Older consumable products listed in this table will continue to be available for purchase until further notice.



SPK 3.0 WGS sample preparation updates



SPK 3.0 WGS sample preparation workflow overview

Comparison of SPK 3.0 to previous SMRTbell library prep Kits

SMRTbell prep kit 3.0 enables simpler, more unified sample prep workflows across different applications

		SMRTbell prep kit 3.0 (NEW)	SMRTbell express template prep kit 2.0
 <p>We recommend using SPK 3.0 for all supported applications</p>			
	Part Numbers	1. SMRTbell prep kit 3.0	1. SMRTbell express TPK 2.0 2. Enzyme cleanup mix 2.0 3. AMPure PB beads
	Samples	24	18
	Time, gDNA to SMRTbell library	4.5 hours	8 hours
	Input DNA / 30-fold human	3 µg	5 µg
	Applications*	WGS, shotgun metagenomics, amplicons, Iso-Seq method, viral sequencing	WGS Amplicons Iso-Seq
	Kit cost	\$1,800	\$1,779
	Sample cost	\$75	\$99

New whole genome sequencing (WGS) library prep protocol using SPK 3.0 enables **reduced workflow times**, uses **reduced DNA input** requirements and **eliminates gel-based size selection** requirements

New HiFi library preparation protocol using SMRTbell prep kit 3.0 for whole genome and metagenomic shotgun sequencing applications

Procedure & checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0 ([102-166-600](https://www.pacb.com/documentation/support/faq)) describes a method for constructing SMRTbell libraries that are suitable for generating HiFi reads on the Sequel II and IISX Systems for **WGS and metagenomic shotgun sequencing applications**.

Procedure Highlights

- Uses **SMRTbell prep kit 3.0** (102-182-70) and supports high-throughput processing using **300 ng – 5 µg** of input genomic DNA amounts
 - We recommend starting with **≥1 µg of input DNA per SMRT Cell 8M** (or ~3 µg for up to a 3 Gb WGS sample to enable running 3 SMRT Cells 8M)
- Multiplexing of samples can be performed using **SMRTbell barcoded adapter plate 3.0** (102-009-200)
- Recommend shearing high-quality gDNA using a **Megaruptor 3 System** (Diagenode)
 - **15 kb – 18 kb** target insert size for large (plant / animal / human) genomes
 - **7 kb – 12 kb** target insert size for small (microbial) genomes
 - **7 kb – 12 kb** target insert size for shotgun metagenomic samples
- **4.5-hour workflow time** to process up to 8 samples from shearing to size selection (6 hours for 24 samples)
 - Time difference is from DNA shearing, which can be performed in sets of 8 samples.
 - Excludes time needed for DNA sizing QC analysis using a Femto Pulse system.
- WGS SMRTbell libraries can be **size-selected using AMPure PB beads** without the need for third-party equipment

Preparing whole genome and metagenome libraries using SMRTbell® prep kit 3.0

PacBio

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing whole-genome sequencing (WGS) libraries from genomic and metagenomic DNA using the SMRTbell prep kit 3.0 for sequencing on PacBio systems.

Overview			
Samples per SMRTbell prep kit 3.0	1–24		
Workflow time	4.5 hours for up to 8 samples; 6 hours for 24 samples Time difference is from DNA shearing, which is done in sets of 8 samples. Excludes measuring DNA size on Femto Pulse system.		
DNA input			
Quantity	300 ng–5 µg per library		
	Human, plant, and animal	Microbes	Metagenomes
DNA size distribution (Femto Pulse system)	50% ≥ 30 kb & 90% ≥ 10 kb	90% ≥ 7 kb	90% ≥ 7 kb
DNA Shearing (Megaruptor 3 system)	Speed 31	Speed 40	Speed 40
Target fragment lengths	15–18 kb	7–12 kb	7–12 kb
Size selection required	AMPure® PB beads	none	none

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
PN 102-166-600 EA V1 18FEB2022

PacBio


PacBio [Documentation](https://www.pacb.com/documentation/support/faq) ([102-166-600](https://www.pacb.com/documentation/support/faq))

APPLICATIONS WHOLE GENOME SEQUENCING

De Novo assembly & variant detection


Microbial assembly

Shotgun metagenomics



SPK 3.0 WGS workflow enables faster library construction times with lower DNA input amounts compared to SMRTbell express TPK 2.0

SPK 3.0 WGS protocol reduces walk-away times and consolidates several enzymatic repair reactions & A-tailing into a single reaction



New SMRTbell prep kit 3.0 consolidates all required core template prep reagents into a single orderable part number


SMRTbell prep kit 3.0 WGS workflow (NEW)

SMRTbell prep kit 3.0
(102-182-700) (NEW)



Included with SPK 3.0 :

- ✓ Library construction reagents (24 Rx)
 - Repair buffer, End repair mix, DNA repair mix
 - SMRTbell adapter
 - Ligation mix, Ligation enhancer
 - Other buffers
- ✓ SMRTbell cleanup beads
- ✓ Nuclease mix



** Note: If performing size selection, use AMPure PB bead size selection kit (102-182-500)


SMRTbell express TPK 2.0 WGS workflow

SMRTbell express TPK 2.0
(100-938-900)



Included with SE TPK 2.0:

- ✓ Library construction reagents (18 Rx)
 - DNA prep buffer, NAD, DNA prep additive, DNA damage repair mix v2, End prep mix
 - SMRTbell overhang adapter v3
 - Ligation mix, Ligation additive, Ligation enhancer
 - Other buffers
- ✗ AMPure PB beads (separate purchase required)
- ✗ SMRTbell enzyme cleanup kit 2.0 (separate purchase required)



New/updated protocols using SPK 3.0 & Binding Kit 3.1/3.2

All WGS applications except ultra-low DNA input now use the same protocol

Protocol consolidation with SPK 3.0 unifies sample prep workflows across different application use cases

Main application	Application subtype / Supported use case	Template prep kit(s)	Binding kit	Procedure & checklist for SMRTbell library prep*	
Whole Genome Sequencing	Large genome WGS, microbial genome WGS, low DNA input & shotgun metagenomics	SPK 3.0	Binding Kit 3.2	Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0 [102-166-600]	1
	Ultra-Low DNA Input Sequencing ¹	SMRTbell express TPK 2.0	Binding Kit 3.2	Preparing HiFi SMRTbell libraries from ultra-low DNA input [101-987-800]	
Viral Sequencing	HiFiViral SARS-CoV-2	SPK 3.0	Binding Kit 3.1	PacBio HiFiViral high-throughput multiplexing for full-viral genome sequencing of SARS-CoV-2 [102-396-100]	2
	AAV Sequencing	SPK 3.0	Binding Kit 3.1	Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 [102-126-400]	3
RNA Sequencing	Iso-Seq Method (Bulk)	SPK 3.0	Binding Kit 3.1 Binding Kit 3.2	Preparing Iso-seq libraries using SMRTbell prep kit 3.0 [102-396-000]	4
	Single-Cell Iso-Seq Method ¹	SMRTbell express TPK 2.0	Binding Kit 3.1 Binding Kit 3.2	Preparing single-cell Iso-seq libraries using SMRTbell express TPK 2.0 [101-892-000]	
Metagenomics	Full-length 16S Sequencing ²	SPK 3.0	Binding Kit 3.1	Preparing multiplexed amplicon libraries using SMRTbell prep kit 3.0 [102-359-000] (using recommended PacBio barcoded 16S primer sequences)	
Targeted Sequencing	Amplicon Sequencing (Barcoded adapters or barcoded gene-specific primers)	SPK 3.0	Binding Kit 3.1 Binding Kit 3.2	Preparing multiplexed amplicon libraries using SMRTbell prep kit 3.0 [102-359-000]	5
	Amplicon Sequencing (Barcoded M13 primers)	SPK 3.0	Binding Kit 3.1 Binding Kit 3.2	Preparing multiplexed amplicon libraries using PacBio barcoded M13 primers and SMRTbell prep kit 3.0 [101-921-300]	6

* Note: SMRTbell express TPK 2.0 protocols are still supported and available:

- Preparing HiFi SMRTbell Libraries using SMRTbell Express Template Prep Kit 2.0 [[101-853-100](#)]
- Preparing Multiplexed Microbial Libraries Using SMRTbell Express Template Prep Kit 2.0 [[101-696-100](#)]
- Preparing HiFi Libraries from Low DNA Input Using SMRTbell Express Template Prep Kit 2.0 [[101-730-400](#)]
- Preparing 10 kb Library Using SMRTbell Express Template Prep Kit 2.0 for Metagenomics Shotgun Sequencing [[101-800-800](#)]
- Preparing SMRTbell Libraries using PacBio Barcoded Universal Primers for Multiplex SMRT Sequencing [[101-791-800](#)]

¹ Ultra-low DNA input and Single-cell Iso-Seq method protocols will be updated to support SPK 3.0 in the future

² Can use multiplexed amplicon SPK 3.0 protocol [102-359-000](#) to prepare **16S** and **HLA** SMRTbell libraries



SPK 3.0 WGS sample preparation workflow details

Procedure & checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0

Procedure & checklist [102-166-600](#) describes a method for constructing SMRTbell libraries using SMRTbell prep kit 3.0 that are suitable for generating high-accuracy long reads on the Sequel II and IIS systems for **whole genome sequencing** (*de novo* assembly, variant detection, microbial multiplexing) and **shotgun metagenomic sequencing**

Protocol Contents

1. Recommendations for **gDNA quantification and sizing QC**.
2. Recommendations for shearing gDNA to the desired target mode size using the **Megaruptor 3** system (Diagenode).
3. Enzymatic steps for preparation of a WGS SMRTbell library using **SMRTbell prep Kit 3.0** (102-182-700). (Instructions for preparing **multiplexed** samples using **SMRTbell barcoded adapter plate 3.0** (102-009-200) are also provided.)
4. Instructions for size-selection of WGS SMRTbell libraries using **AMPure PB bead size selection**. (Size selection is not required for microbial WGS and metagenomic shotgun libraries where retention of shorter fragments is desired.)
5. Guidance for **pooling** barcoded WGS SMRTbell libraries for multiplexed sequencing on a single SMRT Cell.

Preparing whole genome and metagenome libraries using SMRTbell® prep kit 3.0

PacBio

Procedure & checklist

Before you begin

This procedure describes the workflow for constructing whole-genome sequencing (WGS) libraries from genomic and metagenomic DNA using the SMRTbell prep kit 3.0 for sequencing on PacBio systems.

Overview			
Samples per SMRTbell prep kit 3.0	1–24		
Workflow time	4.5 hours for up to 8 samples; 6 hours for 24 samples Time difference is from DNA shearing, which is done in sets of 8 samples. Excludes measuring DNA size on Femto Pulse system.		
DNA input			
Quantity	300 ng–5 µg per library		
	Human, plant, and animal	Microbes	Metagenomes
DNA size distribution (Femto Pulse system)	50% ≥ 30 kb & 90% ≥ 10 kb	90% ≥ 7 kb	90% ≥ 7 kb
DNA Shearing (Megaruptor 3 system)	Speed 31	Speed 40	Speed 40
Target fragment lengths	15–18 kb	7–12 kb	7–12 kb
Size selection required	AMPure® PB beads	none	none

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.
PN 102-166-600 EA V1 18FEB2022

PacBio

PacBio [Documentation](#) (102-166-600)

SPK 3.0 WGS library prep and sequencing workflow timing overview


Go from DNA shearing to sequencing preparation in 1 day with SPK 3.0

Workflow step		Hands-on (min)	Walk-away (min)
DNA shearing*	DNA shearing (1 Megaruptor 3 cycle)	5	45
	1.0X SMRTbell bead cleanup	5	20
SMRTbell library construction*	DNA repair / a-tailing	5	35
	Adapter ligation (barcoded or non-barcoded adapter)	5	30
	1.0X SMRTbell bead cleanup	5	20
	Nuclease treatment	5	15
	AMPure PB bead size selection	10	30
Total		40 min	3.3 hrs

Workflow step		Hands-on (min)	Walk-away (min)
Sequencing preparation (ABC)	Primer annealing (Sequel II primer 3.2)	5	15
	Polymerase binding (Sequel II binding kit 3.2)	5	15
	Complex cleanup (1.2X SMRTbell cleanup beads)	5	20
Total		15 min	0.83 hrs

SPK 3.0 WGS Workflow

(102-166-600)



Input genomic DNA QC recommendations for WGS library construction using SMRTbell prep kit 3.0

- WGS library construction using SMRTbell prep kit 3.0 requires **high-quality, high-molecular weight genomic DNA**.*
- Prior to library preparation, evaluate the **quantity** and **size distribution** of the input gDNA to determine whether it is suitable for the protocol.
- For each input gDNA sample:
 - ❑ Measure concentration and total mass of DNA with a **Qubit High Sensitivity dsDNA Assay** system (Thermo Fisher Scientific)
 - ❑ Measure DNA size distribution with a **Femto Pulse** system (Agilent)
 - ❑ Proceed with SMRTbell library construction if the **gDNA sample quality** is acceptable (see Table below)

Sample type	Input DNA metric	Requirement	Notes
All	Per Library	300 ng – 5 µg	<ul style="list-style-type: none"> • Starting with low DNA input amounts approaching ~300 ng may in some cases produce insufficient amounts of SMRTbell library to load at concentrations that optimize sequencing data yield. • For multiplexing applications, generally aim to use ≥300 ng of DNA input per sample, with a total mass ≥1 µg across all samples
All	Per SMRT Cell 8M	≥1 µg	<ul style="list-style-type: none"> • Start with ≥1 µg of total input DNA per SMRT Cell 8M (for a single sample or across multiple samples when pooling) to enable generation of sufficient library to load at concentrations that optimize sequencing data yield.
Large genome (Animal/plant/human)	Longer than 30 kb	≥50%	<ul style="list-style-type: none"> • Required to achieve target fragment lengths after DNA shearing. • For large genome samples, the Femto Pulse Genome Quality Number (GQN) at 30 kb should be ≥5.0. (Not applicable to microbial and metagenomic samples)
Large genome (Animal/plant/human)	Longer than 10 kb	≥90%	<ul style="list-style-type: none"> • Required for effective AMPure PB bead size selection. • For large genome samples, the GQN at 10 kb should be ≥9.0.
Small genome (microbial/metagenomic)	Longer than 7 kb	≥90%	<ul style="list-style-type: none"> • For microbial and metagenomic samples, the input DNA should be at least as large as the recommended insert lengths of 7–12 kb with a GQN at 7 kb ≥9.0. • Any degradation present should be due to shearing from the extraction process (e.g., bead beating) and not from poor sample handling or storage, or biochemical processes

NEW DNA extraction technical note: Sample preparation for PacBio HiFi sequencing from human whole blood (102-326-500)


Provides best practices for handling human whole blood samples* to generate optimal sequencing performance

Technical note [102-326-500](#) discusses the effect of anticoagulant, sample storage time, storage conditions, and white blood cell count on the sequencing performance of DNA extracted using the Nanobind CBB Big DNA Kit ([NB-900-001-01](#))



Nanobind CBB Big DNA Kit (NB-900-001-01) for isolating HMW DNA from cells, bacteria, & blood.

For optimal HiFi yield and read length performance, store human whole blood samples for fewer than 2 days at 4°C.




Stage	Variable	Best practice for PacBio HiFi sequencing
Before DNA extraction	Sample type	Human whole blood
	Anticoagulant	Potassium EDTA (K ₂ EDTA)
	Sample storage temperature	4 ± 3°C
	Sample storage time	≤ 2 days from collection to extraction
DNA extraction	Volume of whole blood	200 µL
	White blood cell (WBC) count	≥ 4 × 10 ⁶ cells/mL for ≥ 3 µg of DNA
	DNA extraction method	Nanobind CBB Big DNA kit
After DNA extraction	DNA storage	Rest 1 day at ambient temperature, then store at 4 ± 3°C
	DNA size distribution	<ul style="list-style-type: none"> 90% of DNA ≥ 10 kb (genomic quality number at 10 kb ≥ 9.0) 50% of DNA ≥ 30 kb (genomic quality number at 30 kb ≥ 5.0)
	UV absorbance	<ul style="list-style-type: none"> A_{260/280} nm ≥ 1.7 A_{260/230} nm ≥ 1.5

The Femto Pulse System is recommended for DNA sizing QC of genomic DNA for WGS applications

- Femto Pulse System (Agilent) is highly recommended for DNA sizing QC of input genomic DNA and SMRTbell libraries
 - Enables sizing of gDNA samples ranging from 1,300 bp to 165 kb
 - Requires <1 ng of sample DNA
 - Can analyze up to 12 samples in <1.5 hours
- The Femto Pulse system can be used in place of traditional pulse-field gel electrophoresis (PFGE) to quickly assess the initial integrity of genomic DNA, evaluate shears, determine appropriate size-selection thresholds, and conduct final QC before preparing libraries for SMRT Sequencing

Femto Pulse outputs quality metrics such as the **Genomic Quality Number (GQN)*** to quickly score the integrity of HMW gDNA

Ave. Smear Size (bp)	GQN Set at 30 kb
12,147	0
23,339	1.5
45,304	6.4
57,789	7.1
73,267	7.8
94,045	7.8
109,968	8.2
164,292	8.8



Femto Pulse System



Femto Pulse system offers a simplified QC workflow to generate SMRTbell libraries for WGS sequencing in **reduced time**, **and conserves sample** by using femtogram ranges of input DNA

The Megaruptor 3 system is recommended for shearing genomic DNA For WGS applications

- **Megaruptor 3 system** (Diagenode) is highly recommended for DNA shearing*
 - Up to 8 samples can be sheared in parallel in ~45 minutes for high-throughput applications
 - Achieving the same size distribution across multiple samples provides more consistent sequencing performance
- Recommended library insert size distributions and Megaruptor 3 shear speed settings to use for different WGS applications are summarized on Page 7 in the procedure
 - Bring input gDNA to a final volume of **100 – 130 μL with Low TE buffer** [10 mM Tris-HCl (pH 8.0) + 0.1 mM EDTA] to target a DNA concentration of **3 – 39 $\text{ng}/\mu\text{L}$ (ideal: 30 $\text{ng}/\mu\text{L}$)**
 - Perform shearing (**1 cycle**) using the conditions described in the table below

Megaruptor 3 System



Application	Recommended library Insert size (mode)	Recommended Megaruptor 3 shear speed setting
Animal / plant / human WGS	15 kb - 18 kb	31
Microbial WGS or shotgun metagenomics	7 kb - 12 kb	40


Because the response of individual gDNA samples can differ, **optimization of shearing conditions** may be needed to achieve the desired fragment distribution



* **Note:** The g-TUBE (Covaris) device generates a broader DNA fragment size-distribution compared to the Megaruptor 3 system. As a result, HiFi read quality and overall HiFi data yield may be reduced due to the residual presence of very large DNA fragments generated by g-TUBEs. For additional guidance, see [Technical Note: Covaris g-TUBE DNA Shearing for SMRTbell prep kit 3.0 \(102-326-501\)](#) or contact [PacBio Technical Support](#) or your local Field Applications Scientist.

Example Megaruptor 3 shearing results for a human genomic DNA sample


Megaruptor 3 System Speed Setting 31



Femto Pulse DNA sizing QC analysis of a human gDNA sample sheared using a Megaruptor 3 with speed setting 31 (1-cycle shear). The fragment size distribution mode is 17.8 kb.


Example Megaruptor 3 shearing results for microbial genomic DNA samples

Megaruptor 3 System Speed Setting 40



Femto Pulse DNA sizing QC analyses of four different microbial gDNA samples sheared using a Megaruptor 3 System with speed setting 40 (1-cycle shear). The mean sheared DNA fragment size for all samples is ~7 kb – 10 kb.

g-TUBE 3287 x g [7000 RPM with Eppendorf MiniSpin Plus)



Femto Pulse DNA sizing QC analyses of four different microbial gDNA samples sheared using g-TUBES with a centrifugation speed of 3287 x g. The mean sheared DNA fragment size for all samples is ~7 kb – 10 kb.

SMRTbell barcoded adapter plate 3.0 is recommended for barcoding WGS samples

For Sequel II and IIe Systems, SMRTbell barcoded adapter plate 3.0 (102-009-200) is available for multiplexing up to 96 microbes per SMRT Cell 8M.

- Contains 96 barcoded adapters to support multiplexed SMRTbell library construction for up to 96 samples using SPK 3.0
- Can be used for Microbial Assembly and any other WGS or amplicon sequencing application that employs barcoded overhang adapters
- Each barcoded adapter contains a 5 bp padding sequence for more uniform ligation performance across different barcode sequences
- Each well on the plate contains a barcoded adapter with a unique 10-base pair PacBio barcode sequence
- Each barcoded adapter is present in only one well and supports a single reaction
- SMRT Link comes pre-installed with the following barcode set FASTA file containing SMRTbell barcoded adapter plate 3.0 barcode sequences*:
SMRTbell Barcoded Adapter Plate 3.0 (bc2001-bc2096)

	1	2	3	4	5	6	7	8	9	10	11	12
A	BC 2001	BC 2009	BC 2017	BC 2025	BC 2033	BC 2041	BC 2049	BC 2057	BC 2065	BC 2073	BC 2081	BC 2089
B	BC 2002	BC 2010	BC 2018	BC 2026	BC 2034	BC 2042	BC 2050	BC 2058	BC 2066	BC 2074	BC 2082	BC 2090
C	BC 2003	BC 2011	BC 2019	BC 2027	BC 2035	BC 2043	BC 2051	BC 2059	BC 2067	BC 2075	BC 2083	BC 2091
D	BC 2004	BC 2012	BC 2020	BC 2028	BC 2036	BC 2044	BC 2052	BC 2060	BC 2068	BC 2076	BC 2084	BC 2092
E	BC 2005	BC 2013	BC 2021	BC 2029	BC 2037	BC 2045	BC 2053	BC 2061	BC 2069	BC 2077	BC 2085	BC 2093
F	BC 2006	BC 2014	BC 2022	BC 2030	BC 2038	BC 2046	BC 2054	BC 2062	BC 2070	BC 2078	BC 2086	BC 2094
G	BC 2007	BC 2015	BC 2023	BC 2031	BC 2039	BC 2047	BC 2055	BC 2063	BC 2071	BC 2079	BC 2087	BC 2095
H	BC 2008	BC 2016	BC 2024	BC 2032	BC 2040	BC 2048	BC 2056	BC 2064	BC 2072	BC 2080	BC 2088	BC 2096

Figure illustration of mapping between a specific well location and a unique PacBio barcode sequence on a 96-well plate in the SMRTbell barcoded adapter plate ([102-009-200](#))

Reagent kit quantities support a **single use** of each of the 96 barcoded adapters in the plate for SMRTbell library preparations.

Plate Layout (Excel) [[Link](#)]

Barcode Sequences (FASTA) [[Link](#)]



Product insert: SMRTbell barcoded adapter plate 3.0 (96 barcodes, 96 samples) [[Link](#)]



Size-selection of SMRTbell libraries with AMPure PB beads is recommended for WGS applications

- AMPure PB beads are used as the **default size selection method*** to remove short DNA fragments (<5 kb) and enrich for the long fragments when preparing SMRTbell libraries for whole genome sequencing
- AMPure PB bead size selection of SMRTbell templates is performed follows:
 - Prepare a **35% dilution (v/v)** of the AMPure PB bead stock by adding 1.75 mL of resuspended AMPure PB beads to 3.25 mL of Elution Buffer (EB). [35% AMPure PB beads solution can be stored at 4°C for 30 days.]
 - Add 3.1X v/v of resuspended, room-temperature 35% AMPure PB beads solution to each sample tube and incubate for 20 min at RT to allow beads to bind to DNA
 - Place sample tubes on a magnetic rack to immobilize AMPure PB beads; wash samples with 80% ethanol 2X; then elute samples in 15 μ L of EB for 5 min at RT

AMPure PB beads



With high-quality WGS samples, AMPure PB bead size selection can recover sufficient SMRTbell library material to run **up to ~3 or more SMRT Cells 8M per 3 μ g of starting input gDNA**


* **Note:** Although size-selection beads have many advantages, automated DNA size selection systems that utilize gel cassettes offer more flexibility in defining a size cutoff. Three automated DNA size selection tools that may **optionally** be employed for performing size selection on SMRTbell libraries for HiFi WGS applications include the PippinHT, BluePippin and SageELF systems from Sage Science. Note that use of these tools requires **higher input DNA amounts ($\geq 1.5 \mu\text{g}/\text{SMRT Cell 8M}$)**. For more information, refer to [Technical Note: Alternative size selection methods for SMRTbell prep kit 3.0 \(TN103-110921\)](#), which provides detailed guidance for size selection of WGS libraries using automated DNA size selection tools or contact [PacBio Technical Support](#) or your local Field Applications Scientist.

Summary comparison of SPK 3.0 library sample preparation for large genome vs. small genome (microbial / shotgun metagenomic) WGS applications

Large Genome Workflow



Small Genome Workflow



Step	Action or metric	Large genome (Single sample)	Small genome (Multiplexed samples)
1	Input DNA quality	<ul style="list-style-type: none"> 50% ≥ 30 kb & 90% ≥ 10 kb 	<ul style="list-style-type: none"> 90% ≥ 7 kb
	Input DNA amount*	<ul style="list-style-type: none"> Use ≥1 µg of total input DNA per SMRT Cell 8M 	<ul style="list-style-type: none"> Use ≥300 ng of input DNA per sample, with a total mass ≥1 µg across all samples
	Megaruptor system shearing	<ul style="list-style-type: none"> Use speed setting 31 	<ul style="list-style-type: none"> Use Speed setting 40
2	Adapter ligation	<ul style="list-style-type: none"> Use standard (non-barcoded) adapter included with SPK 3.0 	<ul style="list-style-type: none"> Use Barcoded adapter plate 3.0
3	Size selection	<ul style="list-style-type: none"> Use AMPure PB bead size selection 	<ul style="list-style-type: none"> AMPure PB bead size selection is optional; otherwise perform a standard 1X cleanup with SMRTbell cleanup beads
4	Pooling	<ul style="list-style-type: none"> N/A 	<ul style="list-style-type: none"> Perform equal mass pooling with barcoded samples
5	Sequencing prep & data collection	<ul style="list-style-type: none"> 30 hr movie time 	<ul style="list-style-type: none"> 15 hr movie time

* Increase DNA input amounts to ≥1.5 µg per SMRT Cell 8M when using a gel-based automated size selection system.




SPK 3.0 WGS example performance data


Example library QC and sequencing performance for human WGS libraries prepared with SMRTbell prep kit 3.0

SMRTbell library QC and primary sequencing metrics

Size-Selected Library QC




Raw Data Report



A secondary left hand peak may also be visible depending on DNA sample quality

CCS Analysis Report



Input gDNA for Megaruptor 3 shearing	3000 ng
Post-shearing recovery (%)*	2620 ng (87%)
Final yield of AMPure PB bead Size-selected library (%)**	1070 ng (36%)

* Post-shearing recoveries typically ranged from ~70% to >95% when using input human DNA samples (1 µg to 5 µg)

** Final post-size selected library yields typically ranged from ~25% to ~50% when using input human DNA samples (1 µg to 5 µg)

Raw Base Yield	617.65 Gb
Mean Polymerase Read Length	102.8 kb
P0	23.4%
P1	75.0%
P2	1.6%

Example sequencing metrics for a human WGS sample run with Binding Kit 3.2 (Polymerase 2.2) / 85 pM on-plate concentration / 30-h movie time / 2-h Pre-Extension Time / Adaptive Loading Target = 0.85


HiFi Reads	2.7 M
HiFi Base Yield	39.2 Gb
Mean HiFi Read Length	14,490 bp
Median HiFi Read Quality	Q34
HiFi Read Mean # of Passes	12

For SPK 3.0 human WGS libraries, per-SMRT Cell HiFi base yields typically ranged from ~28 to 39 Gb.


Example sequencing performance for human WGS libraries prepared with SMRTbell prep kit 3.0 (cont.)

Primary sequencing metrics (Cont.)

Yield by HiFi Read Length



Read Quality Distribution



HiFi Read Length Summary

Read Length (bp)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ 0	2,707,732	100	39,236,168,651	100
≥ 5,000	2,664,322	98	39,051,919,399	100
≥ 10,000	2,353,137	87	36,541,368,326	93
≥ 15,000	1,164,272	43	21,435,305,025	55
≥ 20,000	294,460	11	6,522,779,501	17
≥ 25,000	21,062	1	559,040,421	1
≥ 30,000	1,012	0	35,294,569	0
≥ 35,000	388	0	15,240,023	0
≥ 40,000	129	0	5,578,841	0

HiFi Read Quality Summary

Read Quality (Phred)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ Q20	2,707,732	100	39,236,168,651	100
≥ Q30	1,811,377	67	25,413,473,886	65
≥ Q40	679,582	25	8,150,599,400	21
≥ Q50	146,257	5	1,355,549,531	3

Example *de novo* assembly performance for human WGS libraries prepared with SMRTbell prep kit 3.0

HiFi WGS data sets generated with SPK 3.0 provide highly contiguous and highly accurate assemblies

HG002 Library ID	Contig_N50_Mbp
64009e_s10_cov30	35.4
64012e_s10_cov30	36.5
64015e_s10_cov30	36.7
64438e_s10_cov30	34.1
64441e_s10_cov30	33.6

HG002 Library ID	deNovo_asm_QV
64009e_s10_cov30	48.3
64012e_s10_cov30	48.3
64015e_s10_cov30	48.2
64438e_s10_cov30	48.3
64441e_s10_cov30	48.2

- Data were generated from five different human HG002 WGS libraries run on five different Sequel IIe systems
- Data were subsampled to 30-fold coverage and assembled using SMRT Link Genome Assembly analysis application

Example variant detection performance for human WGS libraries prepared with SMRTbell prep kit 3.0

HiFi WGS data sets generated with SPK 3.0 provide highly accurate variant calls

HG002 Library ID	INDEL.F1_Score
64009e_s10_cov30	0.995
64012e_s10_cov30	0.994
64015e_s10_cov30	0.994
64438e_s10_cov30	0.993
64441e_s10_cov30	0.994

HG002 Library ID	SNP.F1_Score
64009e_s10_cov30	0.999
64012e_s10_cov30	0.999
64015e_s10_cov30	0.999
64438e_s10_cov30	0.999
64441e_s10_cov30	0.999

- Data were generated from five different human HG002 WGS libraries run on five different Sequel IIe systems
- Data were subsampled to 30-fold coverage and analyzed with [DeepVariant](#)



SMRT Link Sample Setup updates



New High-Throughput Sample Setup mode overview

Specifying multiple samples in one protocol instance

New High-Throughput mode feature provides a more streamlined workflow to efficiently process single samples or multiple samples in parallel using automation

- The user can specify multiple samples to all be processed at once
 - A *single* protocol will be generated for *all* of these samples
 - This greatly *simplifies* batch processing and automation
 - All samples in a batch should have *approximately the same concentration and insert size**
 - This generally requires that the user does a *concentration normalization* step at the end of library prep if working with multiple samples

* **All samples in a group should have substantially equivalent library properties** (i.e., insert sizes and concentrations within +/- 15% of the specified values.)

The user can enter any number of samples (including 1)

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Samples"/>
Comment	<input type="text"/>
Application	HiFi Reads
Binding Kit	Sequel® II Binding Kit 3.2
Number of Samples	<input type="text"/> samples Number of Samples is not entered or is invalid.
SMRT Cells per Sample	<input type="text"/> cells Cells to bind is not entered or invalid.
Available Volume per Sample	<input type="text"/> uL Available volume is not entered or invalid.
Insert Size ⓘ	<input type="text"/> bp Insert size is not entered or invalid.
Sample Concentration ⓘ	<input type="text"/> ng/uL
Cleanup Anticipated Yield ⓘ	<input type="text" value="75"/> %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	<input type="text"/> pM Concentration on plate is not entered or invalid.

Parallel display of multiple batches at sample information entry step

Example: Multiple samples per batch

New Sample Group

	< Sample Group >	< Sample Group >
Actions	Copy Remove Lock Automate	Copy Remove Lock Automate
Name	Batch 1	Batch 2
Comment	Humans 1 2 3 and 4	Mice 1 2 and 3
Application	HiFi Reads	HiFi Reads
Binding Kit	Sequel II Binding Kit 3.2	Sequel II Binding Kit 3.2
Number of Samples	4 samples	3 samples
SMRT Cells per Sample	3 cells	3 cells
Available Volume per Sample	100 uL	100 uL
Insert Size	15000 bp	17500 bp
Sample Concentration	40 ng/uL	60 ng/uL
Cleanup Anticipated Yield	75 %	75 %
Recommended Concentration on Plate	50-90 pM	50-90 pM
Specify Concentration on Plate	50 pM	50 pM
Minimum Pipetting Volume	1 uL	1 uL
Warnings		

All 4 samples in 'Batch 1' group have **equivalent** library properties and are processed with the **same** ABC protocol

All 3 samples in 'Batch 2' group have **equivalent** library properties and are processed with the **same** ABC protocol

Example: Single samples (Each "batch" = 1 sample)

New Sample Group

	< Sample Group >	< Sample Group >	< Sample Group >
Actions	Copy Remove Lock Automate	Copy Remove Lock Automate	Copy Remove Lock Automate
Name	Sample 1	Sample 2	Sample 3
Comment	Humans 1	Human 2	Human 3
Application	HiFi Reads	HiFi Reads	HiFi Reads
Binding Kit	Sequel II Binding Kit 3.2	Sequel II Binding Kit 3.2	Sequel II Binding Kit 3.2
Number of Samples	1 samples	1 samples	1 samples
SMRT Cells per Sample	3 cells	3 cells	3 cells
Available Volume per Sample	100 uL	100 uL	70 uL
Insert Size	20000 bp	17500 bp	16000 bp
Sample Concentration	20 ng/uL	60 ng/uL	20 ng/uL
Cleanup Anticipated Yield	75 %	75 %	75 %
Recommended Concentration on Plate	50-90 pM	50-90 pM	50-90 pM
Specify Concentration on Plate	50 pM	50 pM	50 pM
Minimum Pipetting Volume	1 uL	1 uL	1 uL
Warnings			

Sample 1, Sample 2 and Sample 3 have **non-equivalent** library properties and are processed with three **different** ABC protocols


New Sequel II binding kit component details

Sequel II binding kit 3.1 and cleanup beads (102-333-400)

- Sequel II Primer 3.1 (pre-diluted AND pre-conditioned)
- Sequel II Polymerase 2.1
- SMRTbell Cleanup Beads for complex cleanup
- DNA Internal Control 3.1 (defined 2 kb template bound to Polymerase 2.1)
- Supports ≥ 24 binding reactions, and up to 4 SMRT Cells 8M per binding reaction (96 cells total), depending on use case, sample size and concentration
- Note: Sequel II Binding Kit 3.1 is recommended for inserts < 3 kb

Sequel II Primer 3.1
= Sequencing Primer v4*


DNA Internal Control 3.1
= DNA Internal Control 1.0



Sequel II binding kit 3.2 and cleanup beads (102-333-300)

- Sequel II Primer 3.2 (pre-diluted AND pre-conditioned)
- Sequel II Polymerase 2.2
- SMRTbell Cleanup Beads for complex cleanup
- DNA Internal Control 3.2 (defined 11 kb template bound to Polymerase 2.2)
- Supports ≥ 24 binding reactions, and up to 4 SMRT Cells 8M per binding reaction (96 cells total), depending on use case, sample size and concentration
- Note: Sequel II Binding Kit 3.2 is recommended for inserts > 3 kb

Sequel II Primer 3.2
= Sequencing Primer v5*



Sequel II DNA internal control complex 3.1 and 3.2 descriptions

Sequel II DNA control complexes are fixed SMRTbell templates pre-bound to a polymerase and are used as a spike-in sequencing control


Sequel II DNA Internal Control 3.1



DNA Control 3.1 is the same as DNA Control 1.0

- Included with **Sequel II Binding Kit 3.1**
- Synthetic 2 kb construct bound to Sequel II Polymerase 2.1
- Sequence available in NCBI GenBank (Accession [MG495226](#))

Sequel II DNA Internal Control 3.2



- Included with **Sequel II Binding Kit 3.2**
- Synthetic 11 kb construct bound to Sequel II Polymerase 2.2



- DNA internal control complex is added to the sample at the end of Sample Setup and is intended to serve as an ideal (well-behaved) sample for obtaining good sequencing results on the Sequel IIe System
 - Therefore, poor sequencing results obtained with the control can indicate potential issues with the system, reagents, or consumables
- DNA internal control templates are synthetic constructs with a known sequence and therefore can be easily separated from the sample data
- Note: PacBio **requires** the use of a DNA Internal Control to be eligible for reimbursement requests arising from sequencing run failures

Example Sequel II DNA internal control 3.2 performance

Sequel II DNA Internal Control 3.2 polymerase read length performance more closely matches large insert WGS sample library performance compared to Control 1.0



Sample			Control 1.0				
Well	Name	Movie Time (hrs)	Mean	Poly RL Mean (bp)	Total Reads	Mean	Mode
A01	WGS Sample 1 Pol 2.2	30	113986	56757	6610	0.86	0.89
B01	WGS Sample 2 Pol 2.2	30	108072	55405	6886	0.86	0.89
C01	WGS Sample 3 Pol 2.2	30	111207	52599	4522	0.86	0.89
D01	WGS Sample 4 Pol 2.2	30	110184	55776	6854	0.86	0.89

Sequel II DNA Internal Control 1.0 polymerase read length is typically >40 kb for 30-h movies with WGS samples bound to Polymerase 2.2




Sample			Control 3.2				
Well	Name	Movie Time (hrs)	Mean	Poly RL Mean (bp)	Total Reads	Mean	Mode
A01	WGS Sample 1 Pol 2.2	30	93898	89237	1963	0.88	0.91
B01	WGS Sample 2 Pol 2.2	30	102673	91526	1926	0.89	0.91
C01	WGS Sample 3 Pol 2.2	30	107079	99541	1855	0.89	0.91

Sequel II DNA Internal Control 3.2 polymerase read length is typically >80 kb for 30-h movies with WGS samples bound to Polymerase 2.2



Creating a new sample setup calculation worksheet from the SMRT Link v11.0 Sample Setup homepage

High-Throughput mode is the new default mode to create or edit calculations to define reaction conditions for primer annealing, polymerase binding, complex cleanup and sequencing.



The screenshot shows the PacBio Sample Setup interface. At the top, there is a navigation bar with the PacBio logo and 'Sample Setup' dropdown. Below this, the 'Sample Setup' title is displayed with a gear icon. A row of action buttons includes '+ New Calculation', 'Edit/Print', 'Import', 'Export', and 'Delete'. A 'Version:' dropdown menu is highlighted with a red box, showing 'Classic' and 'High-Throughput' (selected). Below the dropdown, the text 'Sample Setup HT for Sequel II and Sequel IIe' is visible. At the bottom, a table displays a list of calculations with columns for Name, Date Created, Number of Sa..., SMRT Cells per Sa..., Created By, Locked, and Comment. The first row shows 'Restriction Enzyme Di...' with a date of 2022-03-25, 03:05:27 PM, 1 Sa..., and Created By 'ngunathilake'.

High-Throughput mode provides a simplified, streamlined workflow to efficiently process either one sample or multiple samples with similar library properties (such as mean insert size and DNA concentration) in parallel. You can also export the calculated values to a CSV file for laboratory automation.

- HT mode is designed to support higher throughput operation by **simplifying and unifying protocol steps** across multiple samples
- HT mode is designed to be used in an identical manner across **all** applications
- HT mode **only** supports the new Binding Kits 3.1 and 3.2
- HT mode is designed to be automation-compatible
- **Classic mode* does not support Sequel II Binding Kit 3.1/3.2**

High-Throughput Sample Setup mode application type specification

Sample Setup auto-populates application-specific information for selected fields

- Select an **Application Type**
- Once an application is selected, default values are auto-populated for various fields and highlighted in **green**

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	<input type="text" value="Application is not specified."/>
Binding Kit	<input type="text" value="Sequel II Binding Kit 3.1"/> Binding kit is not entered
Number of Samples	<input type="text"/> samples
SMRT Cells per Sample	<input type="text"/> cells
Available Volume per Sample ?	<input type="text"/> uL
Insert Size ?	<input type="text"/> bp
Sample Concentration ?	<input type="text"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="50"/> %
Recommended Concentration on Plate	N/A
Specify Concentration on Plate	<input type="text"/> pM
Minimum Pipetting Volume ?	<input type="text" value="1"/> uL
Warnings	

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	<input type="text" value="Adeno-Associated Virus"/>
Binding Kit	<input type="text"/>
Number of Samples	<input type="text"/>
SMRT Cells per Sample	<input type="text"/>
Available Volume per Sample ?	<input type="text"/>
Insert Size ?	<input type="text"/>
Sample Concentration ?	<input type="text"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="50"/> %
Recommended Concentration on Plate	N/A
Specify Concentration on Plate	<input type="text"/> pM
Minimum Pipetting Volume ?	<input type="text" value="1"/> uL
Warnings	

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	<input type="text" value="HiFi Reads"/>
Binding Kit	<input type="text" value="Sequel II Binding Kit 3.2"/>
Number of Samples	<input type="text"/> samples Number of Samples is not entered or is invalid
SMRT Cells per Sample	<input type="text"/> cells Cells to bind is not entered or invalid.
Available Volume per Sample ?	<input type="text"/> uL Available volume is not entered or invalid.
Insert Size ?	<input type="text"/> bp Insert size is not entered or invalid.
Sample Concentration ?	<input type="text"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="75"/> %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	<input type="text"/> pM Concentration on plate is not entered or invalid.


New AAV application type

High-Throughput Sample Setup mode application type specification (cont.)

Auto-populated fields are highlighted in green color

- The following fields are auto-populated and highlighted in **green**:
 - Binding Kit
 - Cleanup Anticipated Yield

If any auto-populated entry is manually changed to a different value, then the field will be highlighted in **yellow** color



Binding Kit Sequel® II Binding Kit 3.2

Binding Kit Sequel® II Binding Kit 3.1

Sequencing primer version does not need to be specified. Use the sequencing primer included with the Sequel II Binding Kit.

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	My Batch of Samples
Comment	
Application	HiFi Reads
Binding Kit	Sequel® II Binding Kit 3.2
Number of Samples	samples Number of Samples is not entered or is invalid
SMRT Cells per Sample	cells Cells to bind is not entered or invalid.
Available Volume per Sample	uL Available volume is not entered or invalid.
Insert Size	bp Insert size is not entered or invalid.
Sample Concentration	ng/uL
Cleanup Anticipated Yield	75 %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	pM Concentration on plate is not entered or invalid.

Entering sample library information

< Sample Group >	
Actions	<input type="button" value="Copy"/> <input type="button" value="Remove"/> <input type="button" value="Lock"/> <input type="button" value="Automate"/>
Name	My Batch of Samples
Comment	
Application	HiFi Reads
Binding Kit	Sequel II Binding Kit 3.2
Number of Samples	2 samples
SMRT Cells per Sample	3 cells
Available Volume per Sample	15 uL
Insert Size	18000 bp
Sample Concentration	60 ng/uL
Cleanup Anticipated Yield	75 %
Recommended Concentration on Plate	50-90 pM
Warnings	

Note: Using input sample concentrations outside these ranges may lead to lower-than-expected sequencing performance

- Enter the following information for each group of samples:
 - Name**
 - Enter a **sample group name** (if processing multiple similar samples in a group) or a **sample name** (if processing a single sample)
 - Comment**
 - You can **enter information meaningful to you** like a batch identifier for a LIMS and/or information for each sample such as: "1: Sample ID1234", "2: Sample ID5678", etc.
 - Number of Samples**
 - All samples in a group should have **insert sizes and concentrations within +/- 15%***
 - SMRT Cells per Sample**
 - Available Volume per Sample (µL)**
 - Please enter the volume of the **least** abundant sample in this batch, to ensure adequate volume is available for all samples
 - Sample Insert Size (bp)**
 - Enter the mean insert size based on a DNA sizing QC (smear) analysis of the final library. All samples in a batch should be **within +/- 15%*** of the entered Insert Size value
 - Sample Concentration (ng/µl)**
 - All samples in a batch should be **within +/- 15%*** of the entered Sample Concentration value. The acceptable range of input concentrations depends on insert size as shown in the table at right:

Insert size	Concentration range
≥10 kb	20 – 60 ng/µL
3 kb – 9999 bp	6 – 20 ng/µL
1.5 kb – 2999 bp	3 – 10 ng/µL
500 bp – 1499 bp	1 – 3 ng/µL

Entering sample clean-up information

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	<input type="text" value="HiFi Reads"/>
Binding Kit	<input type="text" value="Sequel II Binding Kit 3.2"/>
Number of Samples	<input type="text" value="2"/> samples
SMRT Cells per Sample	<input type="text" value="3"/> cells
Available Volume per Sample ?	<input type="text" value="15"/> uL
Insert Size ?	<input type="text" value="18000"/> bp
Sample Concentration ?	<input type="text" value="60"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="75"/> %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	<input type="text" value="80"/> pM
Minimum Pipetting Volume ?	<input type="text" value="1"/> uL
Warnings	

Sample Setup entry fields

- **Cleanup Anticipated Yield** field is **auto-filled** with default values based on the specified Binding Kit
 - **60%** for Sequel II binding kit 3.1
 - **75%** for Sequel II binding kit 3.2
- Sample cleanup is performed using **SMRTbell cleanup beads** included with the binding kit

Cleanup of the bound complex allows removal of excess polymerase and sequencing primers, resulting in **higher-quality data**

Can adjust **Cleanup Anticipated Yield** based on previous experience.

Sample loading concentration specification

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	HiFi Reads ▼
Binding Kit	Sequel [®] II Binding Kit 3.2 ▾
Number of Samples	<input type="text" value="2"/> samples
SMRT Cells per Sample	<input type="text" value="3"/> cells
Available Volume per Sample ?	<input type="text" value="15"/> uL
Insert Size ?	<input type="text" value="18000"/> bp
Sample Concentration ?	<input type="text" value="60"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="75"/> %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	<input type="text" value="80"/> pM
Minimum Pipetting Volume ?	<input type="text" value="1"/> uL
Warnings	

Sample Setup entry fields

- Specify **Concentration On-Plate (OPLC)* (pM)**

See [Quick reference card – loading and pre-extension recommendations for Sequel II and IIe systems \(101-769-100\)](#) for any specific **changes** to SMRT Link Sample Setup recommendations

Can **adjust OPLC** based on previous experience with a particular sample type.

Minimum pipetting volume field is no longer hidden in Advanced Options

< Sample Group >	
Actions	Copy Remove Lock Automate
Name	<input type="text" value="My Batch of Samples"/>
Comment ?	<input type="text"/>
Application	HiFi Reads ▾
Binding Kit	Sequel II Binding Kit 3.2 ▾
Number of Samples	<input type="text" value="2"/> samples
SMRT Cells per Sample	<input type="text" value="3"/> cells
Available Volume per Sample ?	<input type="text" value="15"/> uL
Insert Size ?	<input type="text" value="18000"/> bp
Sample Concentration ?	<input type="text" value="60"/> ng/uL
Cleanup Anticipated Yield ?	<input type="text" value="75"/> %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	<input type="text" value="80"/> pM
Minimum Pipetting Volume ?	<input type="text" value="1"/> uL
Warnings	

Sample Setup entry fields

- Specify the **Minimum Pipetting Volume**
- This allows you to set a lower limit on pipetting volumes to use in certain protocol steps, such as sample annealing and binding
- We recommend setting this to 1 μL , though in some cases, for example if sample availability is very limited, it may be appropriate to set a value below 1 μL
- Some protocol steps include fixed values of 1 μL that will not be affected by this setting

Note: Sample Setup High-Throughput mode **does not contain an Advanced Options section**

Sample information summary table

Summary table provides complete sample information for review after completion of all required field entries

Sample Group Information

Sample Group	
Name	My Batch of Samples
Comment	
Number of Samples	2
SMRT Cells Per Sample	3
Sample Volume to Use	8.0 uL
Sample Concentration	60 ng/uL 5.13 nM
Insert Size	18000 bp
Binding Kit	Sequel® II Binding Kit 3.2
Cleanup Anticipated Yield	75 %
Concentration On Plate	80 pM
Application	HiFi Reads
Minimum Pipetting Volume	1 uL
Warnings	

- For HT mode, we **no longer prioritize primer:template or polymerase:template ratios**
 - For HT mode, primer concentration is **locked at 20 nM** and the polymerase concentration is **locked at 15 nM**
 - This is more consistent with the biochemistry and allows more efficient use of reagents
- When following the recommended input sample concentration guidance for different insert sizes:
 - Template concentration can range from ~1 nM to 5 nM in the annealing reaction and from ~0.5 nM to 2.5 nM in the binding reaction

Sample Setup HT summary only shows the **starting** sample template concentration and does not display the template concentration in the annealing or binding reactions

Sample setup warnings

Sample Setup flags any incomplete / invalid entries and other sample issues

Available Volume per Sample ⓘ	15 uL
Insert Size ⓘ	bp Insert size is not entered or invalid.
Sample Concentration ⓘ	60 ng/uL
Cleanup Anticipated Yield ⓘ	75 %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	80 pM
Minimum Pipetting Volume ⓘ	1 uL
Warnings	<ul style="list-style-type: none">Insert size is not entered or invalid.

- If any field is **missing** an entry or contains an invalid entry, then **no** calculations will be performed
 - Enter a valid value to enable calculations to be performed

Sample Concentration ⓘ	120 ng/uL Sample concentration should be between 20 and 60 ng/uL
------------------------	---

Sample Concentration ⓘ	5 ng/uL Sample concentration should be between 20 and 60 ng/uL
------------------------	---

- If the sample concentration is **outside** the recommended range, calculations will still be performed **only** if there is a sufficient amount of sample to proceed with annealing and binding for the number of SMRT Cells specified
 - However, annealing and binding reaction conditions may not be optimal and generate lower than expected sequencing performance

Number of Samples	2 samples
SMRT Cells per Sample	9 cells Too many cells requested for the available sample quantity
Available Volume per Sample ⓘ	15 uL Not Enough Sample Volume
Insert Size ⓘ	18000 bp
Sample Concentration ⓘ	60 ng/uL
Cleanup Anticipated Yield ⓘ	75 %
Recommended Concentration on Plate	50-90 pM
Specify Concentration on Plate	80 pM
Minimum Pipetting Volume ⓘ	1 uL
Warnings	<ul style="list-style-type: none">Not Enough Sample VolumeToo many cells requested for the available sample quantity

- If there is an **insufficient** amount of sample to proceed with annealing and binding for the number of SMRT Cells specified, then **no** calculations will be performed
 - Can try reducing the on-plate loading concentration or the number of SMRT Cells to enable calculations to be performed

Importing and exporting sample setup worksheets

Sample Setup supports importing and exporting calculations in CSV format

The screenshot shows the PacBio Sample Setup interface. At the top, there is a header with the PacBio logo, 'Sample Setup', and user information 'smark (Lab Tech)'. Below the header, there is a 'Sample Setup' section with a gear icon. A row of buttons includes '+ New Calculation', 'Edit/Print', 'Import', 'Export', and 'Delete'. The 'Import' and 'Export' buttons are highlighted with a red box. Below the buttons, there is a 'Version:' section with radio buttons for 'Classic' and 'High-Throughput'. The main content area is titled 'Sample Setup HT for Sequel II and Sequel IIe' and contains a table with columns: Name, Date Created, Number of Sa..., and SMRT Cells. A callout box points to the 'Export' button with the following text:



Note: The content of the CSV file generated using the **Export** button in the Sample Setup home screen is **different** from the content of the CSV file generated using the High-Throughput mode's **Automate** button used for lab automation.

	Name	Date Created ↓	Number of Sa...	SMRT Cells
<input checked="" type="checkbox"/>	My Samples	2022-03-25, 09:13:50 PM	2	
<input type="checkbox"/>	Restriction Enzyme ...	2022-03-25, 03:05:27 PM	1	
<input type="checkbox"/>	Restriction Enzyme ...	2022-03-25, 01:43:51 PM	1	ngunathilake false
<input type="checkbox"/>	My Samples	2022-03-25, 10:15:59 AM	10	smark false test
<input type="checkbox"/>	My Samples	2022-03-24, 05:50:35 PM	1	ktran false test

- To **import** a new calculation, first find (or create) a calculation **similar** to that you wish to import, then **export** it in CSV format. You can then customize the exported CSV file as needed, then **import** the modified CSV file

Exporting a sample setup worksheet for lab automation

From the sample setup worksheet, click **Automate** to export the calculated values to a CSV file for lab automation



Final Loading Dilution

For each of your 10 samples, combine the following and protect from light:

Component	High-Throughput Mode Demo	✓	
Prepared Sample	150.0 uL		
Sequel II Loading Buffer 3.2	200.0 uL		
Diluted Internal Control (Dilution 3)	10.0 uL		
Total Volume	360.0 uL		

Load 115 uL of sample per well and/or store at 4C for up to 24 hours before use.*

Row	Variable Name	Example Value
1	Export Version	4
2	Instructions Version	SMRT Link: 11.0.0.145013; Chemistry Bundle: 11.0.0.143406; Params: 11.0.0
3	Sample Name	
4	Request Name	
5	Anneal Temperature (C)	
6	Anneal Time (minutes)	
7	Anneal Volume (uL)	
8	Anneal Dilution	
9	Anneal Buffer Name	
10	Polymerase Name	
11	Sequel II Polymerase Dilution Buffer Volume	171.5
12	Binding Annealed Sample Volume	16
13	Binding Diluted Polymerase Volume	16
14	Binding Incubation Temperature (C)	25
15	Binding Incubation Time (minutes)	15
16	ICD Buffer Name	ABC Buffer
17	ICD1 Buffer Volume	19
18	ICD1 Internal Control Stock Volume	1
19	ICD2 Buffer Volume	19
20	ICD2 Diluted Internal Control (ICD1) Volume	1
21	ICD3 Buffer Volume	95
22	ICD3 Dilution	
23	Cleanup S1 Dilution	
24	Cleanup S2 Dilution	
25	Cleanup S3 Dilution	
26	Cleanup S4 Dilution	
27	Cleanup S5 Dilution	
28	Cleanup S6 Dilution	
29	Cleanup S2 Dilution Total Volume	100
30	Cleanup S3 Bead Solution Volume	120
31	Cleanup S5 Elution Volume	150
32	Final Loading Prepared Sample Volume	150
33	Final Loading Dilution Buffer Volume	200
34	Final Loading Diluted Internal Control (ICD3) Volume	10
35	Final Loading Volume (microliter)	115

The CSV file generated by the **Automate** button in High-Throughput mode includes all the fields that display in the Sample Setup page, with the volumes listed in each table easily accessible for liquid handling automation purposes.

Note: If needed, HT automation scripts will instruct the robot to distribute the final loading dilution volume across the required number of plate wells to prevent liquid overflow.*

SMRT Link Sample Setup ABC (annealing / binding / cleanup) recommendations

Follow SMRT Link Sample Setup instructions for using the recommended reagent kits for each application*

Main application	Application subtype / Supported use case	Template prep kit(s)	Annealing ¹	Binding ^{2,3}	Complex cleanup
Whole genome sequencing	Large genome WGS, microbial genome WGS, low DNA input & shotgun metagenomics	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.2	SMRTbell cleanup beads
	Ultra-low DNA input sequencing	SMRTbell express TPK 2.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.2	SMRTbell cleanup beads
Viral sequencing	HiFiViral SARS-CoV-2	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1	SMRTbell cleanup beads
	AAV sequencing	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1	SMRTbell cleanup beads
RNA sequencing	Iso-Seq method (bulk) ²	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1 Sequel II Binding Kit 3.2	SMRTbell cleanup beads
	Single-cell Iso-Seq method ²	SMRTbell express TPK 2.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1 Sequel II Binding Kit 3.2	SMRTbell cleanup beads
Metagenomics	Full-length 16S sequencing	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1	SMRTbell cleanup beads
Targeted sequencing	Amplicon sequencing (barcoded adapters or barcoded gene-specific primers) ³	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1 Sequel II Binding Kit 3.2	SMRTbell cleanup beads
	Amplicon sequencing (barcoded M13 primers) ³	SPK 3.0	Use sequencing primer included with binding kit	Sequel II Binding Kit 3.1 Sequel II Binding Kit 3.2	SMRTbell cleanup beads

¹ Sequencing primer 3.1 is included with Sequel II binding kit 3.1; Sequencing primer 3.2 is included with Sequel II binding kit 3.2.

² For standard Iso-Seq Method (bulk) and standard Single-Cell Iso-Seq samples, use Sequel II binding kit 3.1. For Iso-Seq Method (bulk) and Single-Cell Iso-Seq samples with a focus on long transcripts (≥3 kb), use Sequel II binding kit 3.2. (**Note:** The default binding kit specified for the Iso-Seq application type in SMRT Link Run Design is Binding kit 3.1.)

³ For amplicons <3 kb, use Sequel II binding kit 3.1. For amplicons ≥3 kb, use Sequel II binding kit 3.2.



**Sample Setup High-Throughput
mode annealing / binding / cleanup
(ABC) and DNA internal control
dilution procedure**

Annealing procedure

Sequencing primer annealing reaction volumes are automatically calculated in SMRT Link Sample Setup

Annealing primer

For each of your 2 samples, combine the following components in a new low-binding tube and pipet to mix.

Component	My Samples	✓
Sample	8.0 uL	
Annealing Buffer	4.0 uL	
Sequel II Primer 3.2	4.0 uL	
Total Volume	16.0 uL	

Annealing Buffer replaces 10x Primer Buffer v2

Incubate at room temperature for 15 minutes then proceed to the next step.

- Sample Setup calculates annealing and binding reaction volumes for each sample after all the entries have been made
- For primer annealing reactions, always **add 1 part buffer and 1 part primer to 2 parts of sample**
- Anneal samples to sequencing primer by incubating at room temperature (15 - 25°C) for 15 min
- After annealing the primer to each sample, proceed to polymerase binding

Use the pre-diluted and pre-conditioned sequencing primer **included** with the Sequel II Binding Kit.* **No primer conditioning step is needed.**

Classic Mode still states: *“Incubate at room temperature for 15 minutes then transfer to a 4C location for immediate use, or store at -20C for long-term use.”*

- Sequel II Binding Kit 3.1 **includes** Sequel II Primer 3.1 (= v4 Primer)
- Sequel II Binding Kit 3.2 **includes** Sequel II Primer 3.2 (= v5 Primer)

Binding procedure

Polymerase binding reaction volumes are automatically calculated in SMRT Link Sample Setup

Polymerase Dilution

Combine the following components in a single low-bind tube and pipet to mix. The prepared volume of diluted polymerase is sufficient to process all specified samples.

Component	My Samples	✓
Polymerase Stock	1.0 uL	
Sequel® II Polymerase Dilution Buffer	49.0 uL	
Total Volume	50.0 uL	

Sequel II Polymerase Dilution Buffer replaces Sequel Binding Buffer

Total volume supports all samples

- All samples require an initial **polymerase dilution step** (i.e., a 50-fold dilution of the polymerase stock tube)
- Diluted polymerase must be used immediately, any remaining should be discarded.
- For binding reactions, always **add 1 part diluted polymerase to 1 part of annealed sample**
- Bind polymerase to annealed samples by incubating at room temperature* for 15 min
- After binding the polymerase to each sample, proceed to complex cleanup

Diluted Polymerase must be used immediately, any remaining should be discarded.

Use the **same tube** from annealing step

We recommend performing a **fresh binding reaction** for each sample whenever possible

Binding

For each of your 2 samples, add Diluted Polymerase and finger tap or pipet to mix.

Component	My Samples	✓
Annealed Sample	16.0 uL	
Diluted Polymerase	16.0 uL	
Total Volume	32.0 uL	

Updated storage time for bound complexes to **4 weeks** (instead of 7 days)

Incubate at room temperature for 15 minutes. Bound complex can be stored at 4C for 4 weeks.(*)

Cleanup procedure

Complex cleanup reaction volumes are automatically calculated in SMRT Link Sample Setup

Purification of Polymerase Bound SMRTbell® Complexes

1. Equilibrate the Clean-up Beads and Sequel II Loading Buffer 3.2 to room temperature.

2. Please add the following buffer volumes to each sample in each batch, as indicated:

	Demo	✓	Notes
# of Samples in Batch	2		
Binding Reaction	32.0 uL		
ABC Buffer	68.0 uL		
Total Volume	100.0 uL		

Sequel II Loading Buffer 3.1 replaces Complex Dilution Buffer (CDB)
Sequel II Loading Buffer 3.2 replaces Adaptive Loading Buffer (ALB)

ABC Buffer replaces Complex Dilution Buffer (CDB)

3. Add the indicated volume of Clean-up Beads to each sample in each batch and gently pipette-mix. Incubate on the benchtop for 10 minutes.

	Demo	✓	Notes
# of Samples in Batch	2		
Clean-up Beads	120.0 uL		

Use 1.2X SMRTbell Clean-up Beads included with the Sequel II Binding Kit

1.2X vol.

4. Place each tube or tube strip in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. Discard the supernatant. DO NOT wash the collected bead pellet with ethanol.

5. Immediately resuspend the beads in the indicated volumes of room temperature Sequel II Loading Buffer 3.2 and pipette-mix:

	Demo	✓	Notes
# of Samples in Batch	2		
Sequel II Loading Buffer 3.2	150.0 uL		


- **Note:** To streamline high-throughput workflows, DNA quantification of pre-cleanup and post-cleanup complexes is **no longer required***
- As a result, the actual sample OPLC in the final dilution step **may vary** from the specified target OPLC
- If needed, users can adjust the Cleanup Anticipated Yield to improve concordance of target OPLC vs. actual OPLC values

• After performing complex cleanup, proceed to prepare the DNA Internal Control

DNA internal control dilution procedure

Sample Setup provides instructions for diluting the DNA Internal Control stock before adding to the final solution for loading into the sample plate well

- Addition of (Spike-In) DNA Internal Control for each sample is **highly recommended**
 - Can help distinguish instrument-related issues from sample library quality issues and enables quicker troubleshooting
- Perform **three** sequential 1:20 dilution steps using ABC Buffer to prepare the DNA Internal Control working solution
- Add 3 μL (per SMRT Cell 8M) of serially diluted DNA Internal Control to final on-plate sample dilution



Note: PacBio **requires** the use of the Sequel II DNA Internal Control for consumables to be eligible for reimbursement requests arising from sequencing run failures

Internal Control Dilution

Perform three sequential dilution steps using ABC Buffer. Use a new low-binding tube for each dilution step.

1. First Dilution

Reagent	Internal Control		
ABC Buffer	19.0 μL		
Sequel II DNA Internal Control Complex 3.2	1.0 μL		

Use the DNA internal control included with the Sequel II Binding Kit

Mix well by flicking tube by hand and pulse spin to collect contents and keep on ice.

2. Second Dilution

Reagent	Internal Control		
ABC Buffer	19.0 μL		
Diluted Internal Control (Dilution 1)	1.0 μL		

Three small-volume dilution steps are performed for better compatibility with automation systems

Mix well by flicking tube by hand and pulse spin to collect contents and keep on ice.

3. Third Dilution

Reagent	Internal Control	✓	
ABC Buffer	19.0 μL		
Diluted Internal Control (Dilution 2)	1.0 μL		

Mix well by flicking tube by hand and pulse spin to collect contents and keep on ice.

Final loading dilution procedure

Follow Sample Setup instructions to add diluted DNA Internal Control and Loading Buffer to your sample

Final Loading Dilution

For each of your 2 samples, combine the following **and protect from light**:

Component	My Samples	✓	
Prepared Sample	100.0 uL		
Sequel II Loading Buffer 3.2	134.0 uL		
Diluted Internal Control (Dilution 3)	6.0 uL		
Total Volume	240.0 uL		

Use the Sequel II Loading Buffer included with the Sequel II Binding Kit

Load 115 uL of sample per well and/or store at 4C for up to 24 hours before use. (*)

(*) Sequencing performance after storage is dependent on DNA sample quality/type and cannot be guaranteed.

Parallel display of multiple batches at annealing and binding steps

Annealing primer

For each sample in each batch, combine the following components in a new low-binding tube (one tube per sample) and pipet to mix.

	Sample 1	Sample 2	My Batch of Samples	✓	
# of Samples in Batch	1	1	1		
Sample	5.0 uL	5.0 uL	12.0 uL		
Annealing Buffer	2.5 uL	2.5 uL	6.0 uL		
Sequel II Primer 3.2	2.5 uL	2.5 uL	6.0 uL		
Total Volume	10.0 uL	10.0 uL	24.0 uL		

Incubate at room temperature for 15 minutes then proceed to the next step.

Total volume of diluted polymerase master mix prepared is sufficient to process **all** samples in **all** batches

Polymerase Dilution

This dilution produces one **Diluted Polymerase Master Mix** with enough volume for all samples in all batches.

Combine the following components in a single low-bind tube and pipet to mix.

Reagent	Diluted Polymerase	✓	
Polymerase Stock	1.0 uL		
Sequel® II Polymerase Dilution Buffer	49.0 uL		
Total Volume	50.0 uL		

Diluted Polymerase must be used immediately, any remaining should be discarded.

Binding

For each sample in each batch, add Diluted Polymerase and finger tap or pipet to mix.

	Sample 1	Sample 2	My Batch of Samples	✓	
# of Samples in Batch	1	1	1		
Annealed Sample	10.0 uL	10.0 uL	24.0 uL		
Diluted Polymerase	10.0 uL	10.0 uL	24.0 uL		
Total Volume	20.0 uL	20.0 uL	48.0 uL		

Incubate at room temperature for 15 minutes. Bound complex can be stored at 4C for 4 weeks.*

Parallel display of multiple batches at cleanup step

Purification of Polymerase Bound SMRTbell® Complexes

1. Equilibrate the Clean-up Beads and Sequel II Loading Buffer 3.2 to room temperature.

2. Please add the following buffer volumes to each sample in each batch, as indicated:

	Humans 1 2 and 3	Humans 4 and 5	Human 6	✓
# of Samples in Batch	3	2	1	
Binding Reaction	20.0 uL	20.0 uL	48.0 uL	
ABC Buffer	80.0 uL	80.0 uL	52.0 uL	
Total Volume	100.0 uL	100.0 uL	100.0 uL	

Cleanup instructions are displayed in table form for **easier viewing**

3. Add the indicated volume of Clean-up Beads to each sample in each batch and gently pipette-mix. Incubate on the benchtop for 10 minutes.

	Humans 1 2 and 3	Humans 4 and 5	Human 6	✓
# of Samples in Batch	3	2	1	
Clean-up Beads	120.0 uL	120.0 uL	120.0 uL	

4. Place each tube or tube strip in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear. Discard the supernatant. DO NOT wash the collected bead pellet with ethanol.

5. Immediately resuspend the beads in the indicated volumes of room temperature Sequel II Loading Buffer 3.2 and pipette-mix:

	Humans 1 2 and 3	Humans 4 and 5	Human 6	✓
# of Samples in Batch	3	2	1	
Sequel II Loading Buffer 3.2	150.0 uL	150.0 uL	150.0 uL	

6. To elute the polymerase-bound complexes, incubate the samples on the benchtop for at least 5 minutes at room temperature.

7. Place each tube or tube strip in a magnetic bead rack until the beads collect to the side of the tube and the solution appears clear.

8. Transfer eluates to new low-binding tubes or tube strips. Place on ice **and protect from light**.

Optional: If desired, you may quantify sample recovery efficiency using a Qubit instrument, but this is not required.

Sample Batches are now Lockable

New Sample Group

	< Sample Group >	< Sample Group >	< Sample Group >
Actions	Copy Remove Lock Automate	Copy Remove Lock Automate	Copy Remove Lock Automate
Name	Humans 1 2 and 3	Humans 4 and 5	Human 6
Comment ?			
Application	HiFi Reads	HiFi Reads	HiFi Reads
Binding Kit	Sequel® II Binding Kit 3.2	Sequel® II Binding Kit 3.2	Sequel® II Binding Kit 3.2

Lock is required before samples can be imported into Run Design module and also sends a finalized version of the instructions to the server for use in Data Set reports.

Sample Setup HT batches can now be imported into Run Design

Run Design / Create New

New Run Design

Run Information

System Type

SEQUEL II SEQUEL IIe

Run Name

Run 04.12.2022 10:03

Sample Information

Import from Sample Setup (Classic only)

Application

Samples

Only locked samples and locked HT groups are available for import.

Search ...

Name	Date Created ↓	Binding Kit	Type	Control Kit	Insert Size	OPLC
<input type="checkbox"/> Humans 1 2 and 3	2022-04-11, 07:57:44 PM	Sequel® II Binding Kit 3.2	HT	Sequel® II D...	12000	50
<input type="checkbox"/> Cypress Test HiFi locked 202...	2022-04-11, 07:11:38 PM	Sequel® II Binding Kit 2.2	Classic	Sequel® II D...	10000	36

Limitation: Import is still performed at the level of **single collections** (individual SMRT Cells) in Run Design.

→ The user will need to import the **same** Sample Setup HT batch worksheet for all pertinent samples and define appropriate sample names (and other run design parameters if needed) for each SMRT Cell to be run



SMRT Link Run Design updates



GUI and new default CCS analysis output changes

Application type specification

Run Design auto-populates application-specific information for selected fields

- If not importing sample information from Sample Setup, start by first selecting an **Application Type**
- Once an application is selected, default values are auto-populated for various fields and highlighted in **green**

Sample Information

▼ SAMPLE 1: , A01, 15 hour movie

Import from Sample Setup

Application Required

Well Sample Name Required

Bio Sample Name Required

Sample Comment

Sample Well

Template Prep Kit Required

Binding Kit Required

Sequencing Kit Required

New AAV application type

Sample Information

▼ SAMPLE 1: , A01, 15 hour movie

Import from Sample Setup

Application Required

Well Sample Name Required

Bio Sample Name Required

Sample Comment

Sample Well

Template Prep Kit Required

Binding Kit Required

Sequencing Kit Required

Application dropdown menu items:

- Whole Genome Sequencing
- HiFi Reads
- Microbial Assembly
- RNA Sequencing
- Iso-Seq Method
- Viral Sequencing
- Adeno-Associated Virus
- Metagenomics
- Full-Length 16S rRNA Sequencing
- Shotgun Metagenomics
- Amplicon Sequencing
- <3kb Amplicons
- >=3kb Amplicons
- Other
- Custom

Sample Information

▼ SAMPLE 1: , A01, 30 hour movie

Import from Sample Setup

Application Required HiFi Reads

Well Sample Name Required

Bio Sample Name Required

Sample Comment

Sample Well A01

Template Prep Kit Required SMRTbell® Prep Kit 3.0

Binding Kit Required Sequel® II Binding Kit 3.2

Sequencing Kit Required Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp)

Application type specification (cont.)

Auto-populated fields are highlighted in green color

- The following fields are auto-populated and highlighted in **green**:
 - Template Prep Kit
 - Binding Kit
 - Sequencing Kit
 - DNA Control Complex
 - Movie Time Per SMRT Cell
 - Pre-Extension Time

Adaptive Loading is automatically enabled when Sequel II binding kit 2.2 or 3.2 is specified as the default binding kit to use for an application

If any auto-populated entry is manually changed to a different value, then the field will be highlighted in **yellow** color

Sample Information

Template Prep Kit <small>Required</small>	SMRTbell® Prep Kit 3.0
Binding Kit <small>Required</small>	Sequel® II Binding Kit 3.2
Sequencing Kit <small>Required</small>	Sequel® II Sequencing Plate 2.0 (4 rxn)
DNA Control Complex	Sequel® II DNA Internal Control Complex 3.2
Insert Size (bp) <small>Required</small>	
Recommended Concentration on Plate (pM)	50-90 pM
On-Plate Loading Concentration (pM) <small>Required</small>	0
Movie Time per SMRT Cell (hours)	30
Use Pre-Extension	<input checked="" type="radio"/> YES <input type="radio"/> NO
Pre-Extension Time (hours)	2
Include 5mC Calls in CpG Motifs	<input checked="" type="radio"/> YES <input type="radio"/> NO

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

Template Prep Kit <small>Required</small>	SMRTbell® Prep Kit 3.0
Binding Kit <small>Required</small>	Sequel® II Binding Kit 3.2



Template Prep Kit <small>Required</small>	SMRTbell® Prep Kit 3.0
Binding Kit <small>Required</small>	Sequel® II Binding Kit 3.1

SMRT Link Run Design setup recommendations

Follow SMRT Link Run Design instructions for using the recommended run conditions for each application

Main application	Application subtype / Supported use case	Template prep kit(s)	Pre-extension time	Movie collection time
Whole genome sequencing	Large genome WGS, Low DNA input WGS & shotgun metagenomics	SPK 3.0	2 hrs	30 hrs
	Microbial WGS	SPK 3.0	2 hrs	15 hrs
	Ultra-low DNA input sequencing	SMRTbell express TPK 2.0	2 hrs	30 hrs
Viral sequencing	HiFiViral SARS-CoV-2	SPK 3.0	0 hrs	8 hrs
	AAV sequencing	SPK 3.0	2 hrs	24 hrs
RNA sequencing	Iso-Seq method (bulk)	SPK 3.0	2 hrs	24 hrs
	Single-cell Iso-Seq method	SMRTbell express TPK 2.0	2 hrs	24 hrs
Metagenomics	Full-length 16S sequencing	SPK 3.0	0.5 hrs	10 hrs
Targeted sequencing	Amplicon sequencing (barcoded adapters or barcoded gene-specific primers)	SPK 3.0	Use default value in Run Design for the specified insert size	10 hrs (<3 kb) 30 hrs (≥3 kb)
	Amplicon sequencing (barcoded M13 primers)	SPK 3.0	Use default value in Run Design for the specified insert size	10 hrs (<3 kb) 30 hrs (≥3 kb)

- Refer to [Quick reference card – Loading and pre-extension time recommendations for the Sequel II and IIE systems \(101-769-100\)](#) for updates to recommended run setup parameters for specific applications

HiFi read generation is automatically enabled in SMRT Link v11.0

IMPORTANT! By default, all newly created Sequel II and IiE system run designs will specify to **automatically perform CCS analysis** and output **only HiFi reads**

Sequel IiE system

- With Sequel IiE system run designs, the **on-instrument CCS (OICCS) analysis workflow** is **automatically** enabled and outputs a `hifi_reads.bam` file that contains **only HiFi (≥QV 20 CCS) reads**.*
 - In the Run Design **Advanced Options** section, users can specify whether to include low quality reads (non-HiFi reads) in the CCS analysis output (i.e., generate the full `reads.bam` file to support legacy CLR experiments). Note that specifying YES for this option **disables** automated on-instrument barcode demultiplexing, 5mC detection and heteroduplex detection, if applicable.

SMRT Link v10.2 Run Design

Pre-Extension Time (hours)

Generate HiFi Reads ON INSTRUMENT IN SMRT LINK
 DO NOT GENERATE

System Type

SEQUEL II SEQUEL IiE

SMRT Link v11.0 Run Design

Pre-Extension Time (hours)

Include 5mC Calls in CpG Motifs YES NO

CCS Analysis will be performed on-instrument to produce HiFi .bam files.

System Type

SEQUEL II SEQUEL IiE

'On-Instrument / In SMRT Link / Do Not Generate' control is now **hidden** and locked to defaults

HiFi read generation is automatically enabled in SMRT Link v11.0 (cont.)

IMPORTANT! By default, all newly created Sequel II and IIE system run designs will specify to **automatically perform CCS analysis** and output **only HiFi reads**

Sequel II system

- With Sequel II system run designs, **CCS analysis in SMRT Link is automatically** enabled and outputs a `hifi_reads.bam` file that contains **only HiFi (\geq QV 20 CCS) reads**.
 - In the Run Design **Advanced Options** section, users can specify whether to include low quality reads (non-HiFi reads) in the CCS analysis output (i.e., generate the full `reads.bam` file to support legacy CLR experiments). Note that specifying YES for this option **disables** automated (in SMRT Link) barcode demultiplexing, 5mC detection and heteroduplex detection, if applicable

SMRT Link v10.2 Run Design

Pre-Extension Time (hours)

Generate HiFi Reads ON INSTRUMENT IN SMRT LINK
 DO NOT GENERATE

System Type

SEQUEL II SEQUEL IIE

SMRT Link v11.0 Run Design

Pre-Extension Time (hours)

Include 5mC Calls in CpG Motifs YES NO

CCS Analysis will be performed in SMRT Link.

System Type


SEQUEL II SEQUEL IIE

'In SMRT Link / Do Not Generate' control is now **hidden** and locked to defaults

HiFi read generation is automatically enabled in SMRT Link v11.0 (cont.)


- The following **HiFi (QV ≥ 20) data-specific** reports are automatically generated in SMRT Link Run QC:
 - Read Length Distribution:** Displays a read length histogram distribution for HiFi reads
 - Read Quality Distribution:** Displays a read quality histogram distribution for HiFi reads
 - Read Length vs Predicted Accuracy:** Displays a heat map of HiFi read lengths and predicted accuracies.

No longer displays 'Other CCS Reads' (<QV 20) or 'Other Reads'




Read Length Distribution
SMRT Link v11.0 Run QC

No longer displays 'Other CCS Reads' (<QV 20)



Raw Quality Distribution
SMRT Link v11.0 Run QC

No longer displays 'Other CCS Reads' (<QV 20)



Predicted Accuracy vs. Read Length
SMRT Link v11.0 Run QC

File and directory structure output by Sequel IIS system ICS v11.0

Example default file and directory structure output by the Sequel IIS system for each SMRT Cell transferred to network storage

```
<your_specified_output_directory>/r64012e_211206ee_183753/1_A01/  
|--m64012e_211206_183753.baz2bam_1.log  
|--m64012e_211206_183753.ccs.log  
|--m64012e_211206_183753.ccs_reports.json  
|--m64012e_211206_183753.ccs_reports.txt  
|--m64012e_211206_183753.consensusreadset.xml  
|--m64012e_211206_183753.hifi_reads.bam  
|--m64012e_211206_183753.hifi_reads.bam.pbi  
|--m64012e_211206_183753.sts.xml  
|--m64012e_211206_183753.zmw_metrics.json.gz  
|--m64012e_211206_183753.transferdone
```

In ICS v11.0, **default** output files include
*.hifi_reads.bam and *.hifi_reads.bam.pbi
instead of *.reads.bam and *.reads.bam.pbi
(Note: Unless low quality reads are optionally included in
CCS output)

If 5mC CpG Detection is performed, the following additional files are output:

```
|-- m64012e_211206_183753.5mc_report.json  
|-- m64012e_211206_183753.primrose.log
```

Advanced option: Include low quality reads in CCS analysis output

Optional: Click **Advanced Options** and specify, for this Run Design **only**, whether to include **low quality reads** (non-HiFi reads) in the CCS analysis output.

- **Default setting = NO** for all application types
- Note that specifying YES for this option **disables** automatic on-instrument barcode demultiplexing, 5mC detection, and heteroduplex insert detection on Sequel IIe systems since lower-quality data are **not** compatible with these automated workflows*

Advanced Options

Use Adaptive Loading YES NO

Immobilization Time (hrs) 2

CCS Analysis Output - Include Low Quality Reads YES NO

CCS Analysis Output - Include Kinetics Information YES NO

Add Data to Project

Maximum Loading Time (hours) 2

Include non-HiFi reads. Doing so will disable automatic demultiplexing, 5mC detection, and heteroduplex insert detection, if applicable.

CCS Analysis Output - Include Low Quality Reads YES NO

CCS Analysis Output - Include Kinetics Information YES NO

Advanced option: Include low quality reads in CCS analysis output (cont.)


Note: Specification in Run Design to include low quality reads (non-HiFi reads) in the CCS analysis output affects on-instrument and in-SMRT Link data processing options

System	Data processing feature	Include low quality reads in CCS analysis output?	
		NO (Default)	YES
Sequel IIe system	CCS analysis	<ul style="list-style-type: none"> • <code>hifi_reads.bam</code> generated on-instrument 	<ul style="list-style-type: none"> • <code>reads.bam</code> generated on-instrument
	Heteroduplex detection	<ul style="list-style-type: none"> • Automated on-instrument: available • In SMRT Link: unavailable 	<ul style="list-style-type: none"> • On-instrument: unavailable • In SMRT Link: unavailable
	5mC detection	<ul style="list-style-type: none"> • Automated on-instrument: available • Automated in SMRT Link: available* 	<ul style="list-style-type: none"> • On-instrument: unavailable • Manual in SMRT Link: available*
	Barcode demultiplexing	<ul style="list-style-type: none"> • Automated on-instrument: available • Automated in SMRT Link: available 	<ul style="list-style-type: none"> • On-instrument unavailable • Manual in SMRT Link: available
Sequel II system	CCS analysis	<ul style="list-style-type: none"> • <code>hifi_reads.bam</code> generated in SMRT Link 	<ul style="list-style-type: none"> • <code>reads.bam</code> generated in SMRT Link
	Heteroduplex detection	<ul style="list-style-type: none"> • Automated in SMRT Link: available 	<ul style="list-style-type: none"> • Manual in SMRT Link: available (<code>subreads.bam</code> input only)
	5mC detection	<ul style="list-style-type: none"> • Automated in SMRT Link: available 	<ul style="list-style-type: none"> • Manual in SMRT Link: available* (<code>subreads.bam</code> or <code>reads.bam</code> input)
	Barcode demultiplexing	<ul style="list-style-type: none"> • Automated in SMRT Link: available 	<ul style="list-style-type: none"> • Manual in SMRT Link: available (<code>reads.bam</code> input only)

Advanced option: Include kinetics information in CCS analysis output

Optional: Click **Advanced Options** and specify, for this Run Design **only**, whether to include kinetics information (used for epigenetics analysis) in the CCS Analysis output.

- **Default setting = NO** for all application types **except** microbial assembly*
- Adjusting this setting in Run Design **overwrites** any previous global setting specified in the SMRT Link configuration home page.
- **Note:** Adding kinetics information can increase the amount of storage used by the output BAM files by up to **5 times**.



Advanced Options

Use Adaptive Loading YES NO

Immobilization Time (hrs)

CCS Analysis Output - Include Low Quality Reads YES NO

CCS Analysis Output – Include Kinetics Information YES NO

Add Data to Project

Note: For **microbial assembly** samples, kinetics information is retained in the CCS analysis output by **default**.




New 5mC detection option

Include 5mC calls in CpG motifs

If selected, kinetic signatures of cytosine bases in CpG motifs will be automatically analyzed to identify the presence of 5mC during on-instrument CCS (Sequel Ii system only) or during CCS analysis in SMRT Link

- **Default setting = YES** when specifying 'HiFi Reads' or 'Custom' application types
- 5mC detection is automatically performed on-instrument with the Sequel Ii system and in SMRT Link with the Sequel II system (data outputs are the **same** for both methods)
- 5mC calls are output in `hifi_reads.bam` as BAM standard MM and ML tags and can be easily visualized in [IGV](#)
- Processing and storage requirements are **minimal**:
 - File size increase is ~5%
 - On-instrument processing time for Sequel Ii systems is ~10 minutes
- Kinetics are not retained in the CCS analysis output by default, but they can **optionally** be retained as well.
- 5mC calls require a **CpG context and symmetrical methylation** (i.e., does not detect hemi-methylated sites)
- Though trained on human data, 5mC detection has been demonstrated to work on non-human data (e.g., plants (Maize)).
- 5mC consensus calling and other tools planned for a future SMRT Link version.
 - For guidance on command line tool options for 5mC analysis, please contact your local PacBio support team or [PacBio Technical Support](#)

Include 5mC Calls in CpG Motifs YES NO





Hypomethylation at transcription start sites

Example IGV plot demonstrating 5mC detection in HiFi reads for a human HG002 sample. Hypomethylation at active transcription start sites can be easily visualized (unpublished data).

Include 5mC calls in CpG motifs (cont.)

- The **5mC CpG detection** utility generates the following reports:
 - **CpG Methylation in Reads:** Plots the cumulative percentage of CpG sites in the sample against the predicted probability of methylation. (Report appears in SMRT Link Run QC and Data Management)
 - **CpG Methylation in Reads (Histogram):** Histogram plot displaying the percentage of CpG sites in the sample versus the predicted probability of methylation (Report appears in Data Management only)






New heteroduplex detection option

Detect and resolve heteroduplex reads

If selected,* any detected heteroduplexes (HD) are separated into separate reads during on-instrument CCS (Sequel IIe system only) or during CCS analysis in SMRT Link






New AAV application type option

Adeno-associated virus (AAV) application type specification

If selected, AAV sequencing mode applies a modified adapter calling algorithm during post-primary analysis and enables automated heteroduplex (HD) read detection and resolution during CCS analysis





New on-instrument demultiplexing option

On-instrument barcode demultiplexing (Sequel IIe system only)

SMRT Link v11.0 Run Design enables specification to perform barcode demultiplexing on-instrument

Demultiplex Barcodes ON INSTRUMENT IN SMRT LINK
 DO NOT GENERATE

- If on-instrument barcode demultiplexing is performed, the following additional files are output:

```
|-- bc1001--bc1001/m64012e_211206_183753.bc1001--bc1001.consensusreadset.xml  
|-- bc1001--bc1001/m64012e_211206_183753.hifi_reads.bc1001--bc1001.bam  
|-- bc1001--bc1001/m64012e_211206_183753.hifi_reads.bc1001--bc1001.bam.pbi  
|-- m64012e_211206_183753.barcodes.fasta  
|-- m64012e_211206_183753.lima.log  
|-- m64012e_211206_183753.lima_counts.txt  
|-- m64012e_211206_183753.lima_guess.json  
|-- m64012e_211206_183753.lima_guess.txt  
|-- m64012e_211206_183753.lima_reports.txt  
|-- m64012e_211206_183753.lima_summary.txt  
|-- m64012e_211206_183753.unbarcoded.consensusreadset.xml  
|-- m64012e_211206_183753.unbarcoded.hifi_reads.bam  
|-- m64012e_211206_183753.unbarcoded.hifi_reads.bam.pbi
```

Each demultiplexed child data set is output into a subfolder containing the same name as the barcode(s), e.g., 'bc1001-bc1001' for easy identification

- **Note:** The **un-demultiplexed hifi_reads.bam file is not** transferred, it is partitioned into the example file structure shown here.
- Command line tool (`lima`) is available to re-merge demultiplexed files to re-create the original un-demultiplexed `hifi_reads.bam` file if needed



SMRT Link Run QC updates



New Instrument Status view

New Instrument Status view

Instrument Status view enables real-monitoring of sequencing ZMW activity across multiple instruments


The screenshot shows the PacBio Run QC interface. At the top, there are tabs for 'RUN STATUS' and 'INSTRUMENT STATUS', with 'INSTRUMENT STATUS' selected. Below this is a table with columns: Instrument Name, Instrument Status, SMRT Cell Status, Run Completion, and Sequencing ZMWs. The table lists five instruments: 64002e (Ready), 64009e (Ready), 64001e (Running), 64263e (Running), and 64012e (Ready with an alarm). To the right of the table, a graph shows the number of active sequencing ZMWs over time for each instrument. A callout box points to the 'INSTRUMENT STATUS' tab, and another callout box explains the data shown in the table.

- From the Run QC home screen, select **Instrument Status**
- Displayed information is refreshed every 30 seconds






Instrument Name ↑	Instrument Status	SMRT Cell Status	Run Completion	Sequencing ZMWs
64002e	Ready	■ ■	Completed 2 hours ago on 2022-03-08, 01:47:44 PM.	3,500,000 3,000,000 2,500,000 2,000,000
64009e	Ready	■ ■ ■ ■	Completed 1 day and 12 hours ago on 2022-03-07, 04:09:07 AM.	3,500,000 3,000,000 2,500,000 2,000,000
64001e	Running	■ ■ ■ ■	In 2 days and 8 hours at 04:13:44 AM.	3,500,000 3,000,000 2,500,000
64263e	Running	■	In 8.4 hours at 12:16:51 AM.	3,200,000 3,000,000 2,800,000 2,600,000
64012e	Ready 🔔	■		

For each instrument connected to the same instance of SMRT Link, this displays the instrument name and its current status, SMRT Cell status, when the run will be completed, any active alarms, and how many sequencing ZMWs are active

Instrument Status and SMRT Cell Status



The screenshot shows the PacBio Run QC interface. At the top, there is a header with the PacBio logo and 'Run QC' with a dropdown arrow. Below the header is a sub-header with a circular icon containing a checkmark and the text 'Run QC'. The main content is a table with three columns: 'Instrument Name ↑', 'Instrument Status', and 'SMRT Cell Status'. The table contains five rows of data.

Instrument Name ↑	Instrument Status	SMRT Cell Status
64002e	Ready	
64009e	Ready	
64001e	Running	
64263e	Running	
64012e	Ready 	

Instrument Status

- A **red alarm symbol** displays next to the instrument status if any errors or warnings appear during a sequencing run.

SMRT Cell Status

- If an instrument does not have a SMRT Cell tray loaded, the SMRT Cell Status field will **not** display any icons.

SMRT Cell Status Icons



Fully green: SMRT Cell has completed sequencing



Half Green: SMRT Cell is in cell prep or SMRT Cell is currently sequencing



White: SMRT Cell is in the queue for sequencing, but cell prep has not started


Run Completion and Sequencing ZMWs

• Run Completion

- Displays the estimated time remaining to complete sequencing run or the time elapsed since the sequencing run completed.
- Also displays the date (in YYYY-MM-DD format) when the last sequencing run was completed
- Note: If an instrument has been idle for 2 days, then the Run Completion field will be blank

• Sequencing ZMWs

- Displays a plot of how many ZMWs on a SMRT Cell are actively sequencing during a movie collection
 - For sequencing runs using **Binding Kit 2.2 and 3.2**, only the number of actively sequencing **singly-loaded ZMWs (P1)** displays
 - For sequencing runs using **Binding Kit 2.1 and 3.1**, the **total number of actively sequencing ZMWs (P1 + P2)** displays
- Note: Both the Y-axis (Sequencing ZMW count) and X-axis (Elapsed Time) are **auto-scaled** as a sequencing run progresses



For a SMRT Cell that achieves $\geq 50\%$ P1 loading and $\geq 10\%$ P0, the ZMW Sequencing plot should typically display a peak value $> 2,000,000$


Note: Elapsed Time axis is **auto-scaled** and so the shape of the plots may change as movie collection progresses

Sequencing ZMWs plot characteristics depend on binding kit type specification

- Sequencing ZMWs (cont.)

- For optimally loaded runs ($P1 > 50\%$) conducted with **Binding Kit 2.2 and 3.2**, the peak concurrent Sequencing ZMWs value will typically be **less than** the final % $P1$ ZMW yield reported in Run QC at the end of a movie collection. (For poorly loaded runs ($P1 < 50\%$), the two values will typically be more similar.)
- For optimally loaded runs ($P1 > 50\%$) conducted with **Binding Kit 2.1 and 3.1**, the peak concurrent Sequencing ZMWs value will typically be **higher** than the final % $P1$ ZMW yield reported in Run QC. (For poorly loaded runs, the two values will typically be more similar.)

Due to terminations, **not all ZMWs are singly-loaded at the same time** – e.g., some ZMWs are singly-loaded only at or near the end of a movie collection, whereas others are singly-loaded only at the beginning (and undergo termination before the end of the movie).



As acquisition progresses, active complexes in **P2** ZMWs will terminate at a certain rate and only when one SMRTbell complex per ZMW is left sequencing will a high-quality (HQ) region within the raw read be detected and the ZMW state will turn to **P1**.

Comparison of Instrument Status view vs. Run Status view

Example instrument status and run status reports for a well-loaded sample achieving >70% *P1* loading

Run QC RUN STATUS INSTRUMENT STATUS

Instrument Name ↑	Instrument Status	SMRT Cell Status	Run Completion	Sequencing ZMWs
64009e	Ready	■ ■ ■ ■	Completed 1 day and 12 hours ago on 2022-03-07, 04:09:07 AM.	

Peak value of Sequencing ZMWs plot for last cell (Sample D01) is ~3,000,000

Run QC RUN STATUS INSTRUMENT STATUS

Expand All


Sample Information >		Run Settings >		Productivity (%)			Reads <										Control >
Well	Name	Movie Time (hrs)	Status	Total Bases (Gb)	P0	P1	P2	HiFi Reads		Polymerase Read Length			Longest Subread		Poly RL Mean (bp)		
								≥Q20 Reads	Yield	Mean Length	Median QV	Mean	N50	Mean	N50		
A01	WGS Sample 1 Pol 2.2	30	Complete	637.19	25.3	72.5	2.2	2817138	36.53 Gb	12967	Q37	109721	224250	17115	21250	98761	
B01	WGS Sample 2 Pol 2.2	30	Complete	531.02	30.5	67.3	2.2	1844770	28.79 Gb	15608	Q32	98521	197750	27164	46250	91774	
C01	WGS Sample 3 Pol 2.2	30	Complete	673.07	15.7	81.8	2.5	2802352	45.45 Gb	16218	Q33	102674	206750	19974	22750	94451	
D01	WGS Sample 4 Pol 2.2	30	Complete	570.69	19.6	76.8	3.6	2178787	31.12 Gb	14282	Q34	92686	190250	24667	40250	87702	

Final *P1* loading for sample D01 is >70% (6,144,000 *P1* ZMWs)

Example Instrument Status report (top) and Run QC report (bottom) for a WGS sample bound with Binding Kit 3.2 and sequenced using a 30-hour movie collection time. The Sequencing ZMWs plot in the Instrument Status report shows that the peak concurrent Sequencing ZMWs value for the final SMRT Cell in the run (Well D01) is approximately 3,000,000 ZMWs, whereas the final %*P1* ZMW yield reported in the corresponding Run QC metrics table for Well D01 is 76.8% (or 6,144,000 *P1* ZMWs).

Comparison of Instrument Status view vs. Run Status view

Example instrument status and run status reports for a poorly-loaded sample achieving <50% P1 loading



Run QC RUN STATUS INSTRUMENT STATUS

Expand All

Sample Information >		Run Settings >		Productivity (%)			Reads <									
Well	Name	Movie Time (hrs)	Status	Total Bases (Gb)	P0	P1	P2	HiFi Reads		Mean Read Length	Longest Subread					
								≥Q20 Reads	Yield	Mean Length	Min QV	Mean	N50	Mean	N50	Poly RL Mean (bp)
A01	WGS Sample 1 Pol 2.2	30	Complete	605.22	26.6	71.6	1.8	2914608	30.72 Gb	10479	Q39	105...	222250	14637	14750	86336
B01	WGS Sample 2 Pol 2.2	30	Complete	572.02	29.8	68.3	1.9	2535778	26.57 Gb	10479	Q39	104...	214750	20328	42750	87317
C01	WGS Sample 3 Pol 2.2	30	Complete	68.73	84.8	14.7	0.5	575681	1.61 Gb	2797	Q30	60179	229750	10683	67750	75807
D01	WGS Sample 4 Pol 2.2	10	Complete	69.41	73.4	25.8	0.8	1199084	2.65 Gb	2210	Q52	34005	61750	5368	17250	29914

Final P1 loading for sample D01 is <30% (2,065,000 P1 ZMWs)


Example Instrument Status report (top) and Run QC report (bottom) for a WGS sample bound with Binding Kit 3.2 and sequenced using a 10-hour movie collection time. The Sequencing ZMWs plot in the Instrument Status report shows that the peak concurrent Sequencing ZMWs value for the final SMRT Cell in the run (Well D01) is approximately 2,000,000 ZMWs, whereas the final %P1 ZMW yield reported in the corresponding Run QC metrics table for Well D01 is 25.8% (or 2,065,000 P1 ZMWs).



New 5mC-specific report plots

New Run QC report plot: 5mC Detections

Plots



- **CpG Methylation in Reads:** Plots the cumulative percentage of CpG sites in the sample against the predicted probability of methylation. (Report appears in SMRT Link Run QC and Data Management)



SMRT Analysis updates



Updated GUI nomenclature

New conceptual structure in SMRT Analysis v11.0 for specifying workflow type

Data processing workflows are now separated into 'Analyses' and 'Utilities'

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Job Name Required
SMRT Analysis Demo - Creating a New Analysis Job

Workflow Type
 ANALYSIS AUTO ANALYSIS DATA UTILITY

Analysis of Multiple Data Sets
One Analysis per Data Set - Identical Parameters

Choose an option when multiples Data Sets are selected.

Datasets

Displaying rows 1 to 11 out of 1364

	Data Set Details >			Sample Details				
<input type="checkbox"/>	Name	Demultiplexed Subsets	Date Created	Well Sample Name	Run Name	Created By	Bio Sample Name	Barcode Name
<input type="checkbox"/>	Rhino_Verif_HG002_...		2022-03-27, 10:59:...	unknown	unknown	obanerjee	HG002	

- **Analysis:** An analysis uses applications designed to produce biologically-meaningful results. These analysis applications **only** accept HiFi reads
- **Data Utilities:** Data processing utilities are used as intermediate steps to producing biologically-meaningful results. Some data utilities accept **only** HiFi reads whereas other data utilities accept **only** subreads (formerly known as “Continuous Long Reads” in previous SMRT Link versions)

Analysis Applications

SMRT Link Analysis Applications are designed to produce biologically-meaningful results

SMRT Link analysis applications accept only HiFi reads as input.

- **Genome Assembly**

- Generate de novo assemblies of genomes, using HiFi reads.

- **HiFi Mapping**

- Align (or map) reads to a user-provided reference sequence.

- **HiFiViral SARS-CoV-2 Analysis**

- Analyze multiplexed viral surveillance samples for SARS-CoV-2, using HiFi reads.

- **Iso-Seq Analysis**

- Characterize full-length transcript isoforms, using HiFi reads.

- **Microbial Genome Analysis**


- Note: This combines and replaces the Microbial Assembly and Base Modification Analysis applications in the previous release.
- Generate de novo assemblies of small prokaryotic genomes between 1.9-10 Mb and companion plasmids between 2 – 220 kb, and identify methylated bases and associated nucleotide motifs.
- Optionally include identification of 6mA and 4mC modified bases and associated DNA sequence motifs.

- **Minor Variants Analysis**

- Identify and phase minor single nucleotide substitution variants in complex populations.

- **Structural Variant Calling**

- Identify structural variants (Default: ≥ 20 bp) in a sample or set of samples relative to a reference.



Data Utilities

PacBio Data Utilities are used as intermediate steps to producing biologically-meaningful results

The following data utilities accept only HiFi reads as input:

- **5mC CpG Detection**
 - Analyze the kinetic signatures of cytosine bases in CpG motifs to identify the presence of 5mC. (Sequel II only.)
- **Demultiplex Barcodes**
 - Separate reads by barcode.
- **Export Reads**
 - Export HiFi reads that pass filtering criteria as FASTA, FASTQ and BAM files.
 - For barcoded runs, you must first run the Demultiplex Barcodes application to create BAM files before using this application.
- **Mark PCR Duplicates**
 - Remove duplicate reads from a HiFi reads Data Set created using an ultra-low DNA sequencing protocol.
- **Trim Ultra-Low Adapters**
 - Trim PCR Adapters from a HiFi reads Data Set created using an ultra- low DNA sequencing library.

The following data utilities accept only subreads as input:

- **Circular Consensus Sequencing (CCS)**
 - Identify consensus sequences for single molecules.





Updated Microbial Genome Analysis application

Updated Microbial Genome Analysis application

Use SMRT Link v11.0 Microbial Genome Analysis application to perform microbial assembly and base modification detection using HiFi reads

The screenshot shows the PacBio SMRT Analysis interface. At the top, it says 'PacBio SMRT Analysis' and 'Create New Analysis'. Below that, there are two steps: '1. Select Data' and '2. Select Analysis'. The '2. Select Analysis' step is active. A dropdown menu is open, showing 'Microbial Genome Analysis' selected. Below the dropdown, there are buttons for 'Import Analysis Settings' and 'Export Analysis Settings'. Further down, there are two sections: 'Run Base Modification Analysis' and 'Find Modified Base Motifs'. Both sections have radio buttons for 'ON' and 'OFF', with 'ON' selected. At the bottom, there is a button for 'Advanced Parameters'.


- Generate *de novo* assemblies of small prokaryotic genomes between 1.9-10 Mb and companion plasmids between 2 – 220 kb, and identify methylated bases and associated nucleotide motifs.
- Optionally include identification of 6mA and 4mC modified bases and associated DNA sequence motifs. (This requires kinetic information.)
 - Unlike 5mC calling, microbial base modification detection is performed off-instrument (i.e., in SMRT Link only)
 - This requires a Run Design to specify that kinetic information be retained in the CCS analysis output
 - For the Microbial Assembly application type, Run Design automatically defaults to specifying YES for the 'CCS Analysis Output – Include Kinetics Information' field

The screenshot shows a dropdown menu for 'Application' with 'Microbial Assembly' selected. Below the dropdown, there is a radio button for 'CCS Analysis Output – Include Kinetics Information' set to 'YES'.

Note: This combines and **replaces** the Microbial Assembly and Base Modification Analysis applications in the previous SMRT Link release.

Updated Microbial Genome Analysis application (cont.)

View microbial assembly results, detected base modifications and identified modified base motifs in a single analysis job report





APPENDIX: Technical documentation & applications support resources

HiFi sequencing and software v11.0 release technical documentation

Sequel I/II system documentation

- Sequel II and Sequel I/II systems operations guide ([101-774-700](#))
- Sequel II/I/II system v11.0 release notes ([102-279-600](#))
- Sequel II and I/II systems: Data files ([102-144-100](#))
- Sequel I/II system: Location of HiFi reads files ([102-110-200](#))
- Quick reference card – Loading and pre-extension recommendations for the Sequel II and I/II systems ([101-769-100](#))
- Pacific Biosciences glossary of terms ([000-710-267](#))

SMRT Link & other data analysis documentation

- Brief primer and lexicon for PacBio SMRT sequencing webpage ([v11.0](#))
- PacBio bioinformatics file formats documentation webpage ([v11.0](#))
- SMRT Link v11.0 cloud reference guide ([102-295-600](#))
- SMRT Link v11.0 release notes ([102-279-500](#))
- SMRT Link v11.0 software installation guide ([102-278-600](#))
- SMRT Link v11.0 user guide ([102-278-200](#))
- SMRT Link v11.0 web services API use cases ([102-298-700](#))
- SMRT Tools v11.0 reference guide ([102-278-500](#))

HiFi sequencing and software v11.0 release technical documentation (cont.)

Application technical overviews

- HiFi sequencing and software v11.0 release: Technical overview for Sequel II & Sequel IIe system users ([102-399-900](#))
- Technical overview: Adeno-associated virus (AAV) library preparation using SMRTbell prep kit 3.0 ([102-390-400](#))
- Technical overview: Iso-Seq library preparation using SMRTbell prep kit 3.0 ([102-393-400](#))
- Technical overview: Multiplexed amplicon library preparation using SMRTbell prep kit 3.0 ([102-395-900](#))
- Technical overview: Multiplexed SARS-CoV-2 library preparation for full-viral genome sequencing using SMRTbell prep kit 3.0 ([102-399-300](#))
- Technical overview: Whole genome and metagenome library preparation using SMRTbell prep kit 3.0 ([102-390-900](#))

DNA extraction literature

- Circulomics Nanobind [Application notes](#), [Kit handbooks](#), and [Protocols](#)
- Technical note: Preparing DNA for PacBio HiFi sequencing – Extraction and quality control ([102-193-651](#))
- Technical note: Preparing samples for PacBio whole genome sequencing for *de novo* assembly – collection and storage ([TN100-040518](#))
- Technical note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))

HiFi sequencing and software v11.0 release technical documentation (cont.)

Sample preparation literature

- Overview – Sequel systems application options and sequencing recommendations ([101-851-300](#))
- Procedure & checklist – Preparing multiplexed AAV SMRTbell libraries using SMRTbell prep kit 3.0 ([102-126-400](#))
- Procedure & checklist – Preparing multiplexed amplicon libraries using SMRTbell prep kit 3.0 ([102-359-000](#))
- Procedure & checklist – Preparing multiplexed amplicon libraries using PacBio barcoded M13 primers and SMRTbell prep kit 3.0 ([101-921-300](#))
- Procedure & checklist – PacBio HiFiViral high-throughput multiplexing for full-viral genome sequencing of SARS-CoV-2 using SMRTbell prep kit 3.0 ([102-396-100](#))
- Procedure & checklist – Preparing Iso-Seq libraries using SMRTbell prep kit 3.0 ([102-396-000](#))
- Procedure & checklist – Preparing whole genome and metagenome sequencing libraries using SMRTbell prep kit 3.0 ([102-166-600](#))
- Quick reference card – Loading and pre-extension recommendations for the Sequel II and Ile systems ([101-769-100](#))
- Technical note: Alternative size selection methods for SMRTbell prep kit 3.0 ([TN103-110921](#))
- Technical note: Covaris g-TUBE DNA shearing for SMRTbell prep kit 3.0 ([102-326-501](#))

Example PacBio data sets

5mC detection at CpG sites	Dataset	Data type	PacBio system
5mC detection at CpG sites	Human HG002 CpG methylation status	HiFi Reads	Sequel Ile System



www.pacb.com

Research use only. Not for use in diagnostic procedures. © 2022, PacBio. All rights reserved. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omnione, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of Pacific Biosciences of California, Inc. (PacBio). All other trademarks are the sole property of their respective owners.